

## Anodisc Membrane Filter

## 알루미늄 옥사이드, 멤브레인 필터, 아노디스크

Cat. No.	Model	Diameter	Pore Size	PK Price ₩
F13-07-475	68097003	Ø13 mm	0.02 µm	100ea 1,074,000
F13-07-484	68097013		0.10 µm	100ea 942,000
F13-07-493	68097023		0.20 µm	100ea 975,000
F13-07-502	68096002	Ø25 mm	0.02 µm	50ea 916,000
F13-07-511	68096012		0.10 µm	50ea 896,000
F13-07-520	68096022		0.20 µm	50ea 839,000
F13-07-529	68095002	Ø47 mm	0.02 µm	50ea 1,570,000
F13-07-538	68095012		0.10 µm	50ea 1,869,000
F13-07-547	68095022		0.20 µm	50ea 1,589,000

- 친수성 (Hydrophilic)
- Compatible with most solvent and aqueous material
- High pore density and narrow pore size distribution
- Whatman



## Silver Membrane Filter

## 실버 멤브레인 필터

Cat. No.	Model	Diameter	Pore Size	PK Price ₩
F13-157-096	45330	Ø13 mm	0.20 µm	100ea 1,331,100
F13-157-104	45329		0.45 µm	100ea 1,331,100
F13-157-112	45328		0.80 µm	100ea 1,331,100
F13-157-120	45327		1.20 µm	100ea 1,331,100
F13-157-128	45326		3.00 µm	100ea 1,331,100
F13-157-136	45325		5.00 µm	100ea 1,331,100
F13-157-144	45336	Ø25 mm	0.20 µm	50ea 778,800
F13-157-152	45335		0.45 µm	50ea 778,800
F13-157-160	45334		0.80 µm	50ea 778,800
F13-157-168	45333		1.20 µm	50ea 778,800
F13-157-176	45332		3.00 µm	50ea 778,800
F13-157-184	45331		5.00 µm	50ea 778,800
F13-157-192	45348	Ø47 mm	0.20 µm	25ea 984,400
F13-157-200	45347		0.45 µm	25ea 984,400
F13-157-208	45346		0.80 µm	25ea 984,400
F13-157-216	45345		1.20 µm	25ea 984,400
F13-157-224	45344		3.00 µm	25ea 984,400
F13-157-232	45343		5.00 µm	25ea 984,400

- 내화학성 우수, 427°C 내열, 재사용 가능
- Molecular Bonded Filter Media, 비흡수성 및 흡수성
- Used as the Collection Media and Subsequent X-ray Diffraction Substrate for Quantifying Unknown Materials and Compounds
- Sterlitech
- 재질 : 99.07% Pure Silver



## PTFE Filtering Membrane

## 테프론 필터 페이퍼

Cat. No.	Model	지름	Pore Size	PK Price ₩
F12-460-081	N1690-08	Ø 13mm	0.05 µm	10ea 54,000
F12-460-085	N1690-28	Ø 25mm	0.05 µm	10ea 71,900
F12-460-089	N1690-48	Ø 47mm	0.05 µm	10ea 117,300
F12-460-093	N1690-52		0.2 µm	10ea 115,100
F12-460-097	N1690-56		0.45 µm	10ea 135,900
F12-460-101	N1690-60		1.0 µm	10ea 116,600
F12-460-105	N1690-64		5.0 µm	10ea 116,300
F12-460-109	N1690-80		1.0 µm	10ea 486,100
F12-460-113	N1690-84	Ø 90mm	5.0 µm	10ea 484,100

- 재질 : PTFE , -200°C~250°C 내열
- Thickness : 0.2mm
- BOLA

