


PTFE Filtering Tile
PTFE 테프론 필터 플레이트

Cat. No..	Model	Pore Size	Thickness	EA Price ₩
F12-460-011	N1610-10	5 μm	1mm	1,086,500
F12-460-015	N1616-10	10 μm	1mm	1,086,500
F12-460-019	N1616-20	10 μm	2mm	1,453,600
F12-460-023	N1616-30	10 μm	3mm	3,264,000



- 재질 : PTFE , MicroPorous PTFE , -200 °C to 250°C 내열
- 테프론 재질로 최상의 내화학성
- Size : 320 x 320mm
- BOLA

Membrane Filter for RO, NF, UF
RO, NF, UF 용 멤브레인 필터, Ø47 mm

- Pore Size
 - RO (Reverse Osmosis) : 0.001 ~ 0.0001 um (해수의 담수화)
 - NF (Nano Filtration) : 0.005 ~ 0.001 um
 - UF (Ultra Filtration) : 0.01 ~ 0.001 um
- RO, NF, UF용 High Pressure Filter Holder는 437 Page 참조
- Sterlitech
- 분리 대상
 - RO (Reverse Osmosis) : 350 MWCO 이하
 - NF (Nano Filtration) : 350 ~ 1,000 MWCO
 - UF (Ultra Filtration) : 1,000 ~ 300,000 MWCO

Reverse Osmosis / RO 멤브레인

Cat. No.	Model	제조사	Polymer	Pore Size	Rej-size	pH Range	Typical Flux / psi GFD	PK Price ₩
F13-157-248	AG-475	GE	Polyamide	0 MWCO	99.5	4 ~ 11	26 / 225	5ea 171,600
F13-157-256	AK-475				99		26 / 115	5ea 171,600
F13-157-278	SE-475		TF (Thin Film)		98.9	2 ~ 11	22 / 425	5ea 171,600

Nano-filtration / NF 멤브레인

Cat. No.	Model	제조사	Polymer	Pore Size	Rej-size	pH Range	Typical Flux / psi GFD	PK Price ₩
F13-157-320	DK-475	GE	TF (Thin Film)	0 MWCO	98-MgSO ₄	2 ~ 11	22 / 100	5ea 171,600
F13-157-328	DL-475				96-MgSO ₄		31 / 105	5ea 171,600

Ultra-filtration / UF 멤브레인

Cat. No.	Model	제조사	Polymer	Pore Size	Rej-size	pH Range	Typical Flux / psi GFD	PK Price ₩
F13-157-392	GE-475	GE	Polyamide	-	1000 Dalton	2 ~ 11	18 / 400	5ea 171,600
F13-157-400	GH-475		TF (Thin Film)	1000 MWCO	1 K-PEG		20 / 150	5ea 171,600
F13-157-408	GK-475			2000 MWCO	2 K-PEG		17 / 75	5ea 171,600
F13-157-440	PT-475		Polyethersulfone	5000 MWCO	5 K	-	90 / 50	5ea 171,600
F13-157-448	PW-475			20000 MWCO	20 K	2 ~ 11	85 / 30	5ea 171,600
F13-157-464	MW-475			Ultrafilic	100000 MWCO	100 K	1 ~ 10	176 / 20

PTFE Filtering Rod
PTFE 테프론 필터 로드

Cat. No..	Model	Pore Size	Size	EA Price ₩
F15-460-011	N1505-28	5 μm	od 28 × L100 mm	158,000
F15-460-015	N1505-32	5 μm	od 32 × L120 mm	205,200
F15-460-019	N1510-28	10 μm	od 28 × L100 mm	198,000
F15-460-023	N1520-28	50 μm	od 28 × L100 mm	151,000



- 재질 : MicroPorous PTFE , -200 °C to 250°C 내열
- BOLA