

## Thimble Filter 원통 여과지

### • Whatman

Grade	Pore Size	Thickness	재질	Max. Temp.
603	10 um	1 mm	Cellulose	120°C
603G	0.8 um	1.5 ~ 2 mm	Glass Fiber	500°C
603Q	-	2 mm	Quartz Fiber	1000°C

### Grade 603 : Pore 10 um, 두께 1 mm

Cat. No.	Model	Size (od×id)	Height	PK Price ₩
F15-07-203	2800266	28×26 mm	60 mm	25ea 183,000
F15-07-205	2800288	30×28 mm	80 mm	25ea 302,000
F15-07-209	2800280	30×28 mm	100 mm	25ea 296,000
F15-07-215	2800308	32×30 mm	80 mm	25ea 302,000
F15-07-221	2800338	35×33 mm	80 mm	25ea 191,000
F15-07-229	2800373	39×37 mm	130 mm	25ea 408,000
F15-07-237	2800432	45×43 mm	123 mm	25ea 234,000
F15-07-244	2800608	62×60 mm	180 mm	25ea 773,000

### Grade 603G : Pore 0.8 um, 두께 1.7 mm

Cat. No.	Model	Size (od×id)	Height	PK Price ₩
F15-07-363	2814259	29×25 mm	90 mm	25ea 504,000
F15-07-372	2814300	33×30 mm	100 mm	25ea 471,000
F15-07-379	2814432	46×43 mm	123 mm	25ea 494,000

### Grade 603Q : 두께 2 mm

Cat. No.	Model	Size (od×id)	Height	PK Price ₩
F15-07-410	2812259	29×25 mm	90 mm	10ea 169,000



## Thimble Filter 원통 여과지

### • Advantec

Grade	Pore Size	Thickness	재질	Max. Temp.
No.84	8 um	1.5 mm	Cellulose	120°C
No.86R	1 um	1.6 mm	Glass Fiber	500°C
No.88R	-	2.2 mm	Quartz Fiber	1000°C
No.88RH	-	2.2 mm		

- No. 88RH : Dust 농도 측정계량시 감량을 줄이기 위해 900°C 에서 열처리한 여과지

### Grade 84 : Pore 8 um, 두께 1.5 mm (Cellulose)

Cat. No.	Model	Size (od×id)	Height	PK Price ₩
F15-510-333	35400281	28×25 mm	100 mm	25ea 69,900
F15-510-337	35400331	33×31 mm	120 mm	25ea 77,200

### Grade 86R : Pore 1 um, 두께 1.6 mm (Glass Fiber)

Cat. No.	Model	Size (od×id)	Height	PK Price ₩
F15-510-340	35600280	28×24 mm	100 mm	25ea 144,600
F15-510-344	35600330	33×29 mm	120 mm	25ea 169,900

### Grade 88R : 두께 2.2 mm (Quartz Fiber)

Cat. No.	Model	Size (od×id)	Height	PK Price ₩
F15-510-347	35810250	25×22 mm	90 mm	10ea 83,200

### Grade 88RH : 두께 2.2 mm (900 열처리, Quartz Fiber)

Cat. No.	Model	Size (od×id)	Height	PK Price ₩
F15-510-351	35820251	25×22 mm	90 mm	10ea 92,800

