


Filter Holder, Syringe Filter Type 주사기 타입 필터 홀더, SUS

Cat. No.	Model	Description	Application	EA Price ₩
F03-07-101	301000	up to 10 ml	13 mm Membrane 용	309,600
F03-07-108	301200	up to 50 ml	25 mm Membrane 용	527,200

- 재질 : Stainless Steel 304, PTFE O-ring & Gasket
- Syringe Filter 형식으로 주사기와 같이 사용하기에 적합
- Female Luer Inlet, Male Luer Nozzle Outlet • Re-usable
- Sterlitech



PSF Filter Holder, Syringe Filter type 주사기 타입 필터 홀더, PSF

Cat. No.	Model	재질	Application	EA Price ₩
F03-46-101	T521090	PSF	25 mm	48,700

- 재질 : PSF, -10 ~ 165°C 내열, Autoclavable
- Syringe Filter 형식으로 주사기와 같이 사용하기에 적합
- Re-usable, 반복 멸균하여 사용 가능
- Female Luer Inlet, Male Luer Nozzle Outlet • Tarsons



In-line Filter Holder 인 라인 필터 홀더

- Sterlitech

PP In-line and aerosol Holder

Cat. No.	Model	Description	Application	EA Price ₩
F03-157-129	43303010	PP Filter Holder	25 mm	28,300
F03-157-136	43303020		47 mm	107,000

- PP Body and Silicon O-ring



PP

PFA In-line Filter Holder

Cat. No.	Model	Description	Application	EA Price ₩
F03-157-143	43307010	PFA Filter Holder	25 mm	270,200
F03-157-150	43307020		47 mm	628,900

- PFA Body and O-ring



PFA

Stainless Steel In-line Filter Holder

Cat. No.	Model	Description	Application	EA Price ₩
F03-157-157	357000	Stainless Steel Filter Holder	25 mm	949,900
F03-157-164	17307200		47 mm	1,826,700

- Stainless Steel 304SS Body, PTFE Gasket and Silicon O-ring



Stainless Steel

Stainless Steel In-line Filter Holder, Vented

Cat. No.	Model	Description	Application	EA Price ₩
F03-157-171	357400	SUS 304, Vented	47 mm	2,472,000

- Stainless Steel 304SS, PTFE Gasket and FPM O-ring



Stainless Steel-Vented

In-line Filter Holder, PSF PSF 인 라인 필터 홀더

Cat. No.	Model	재질	Application	EA Price ₩
F03-46-111	T521080	PSF	47 mm	80,500

- 재질 : PSF, -10 ~ 165°C 내열, Autoclavable
- Liquid나 Gas의 Particle이나 Microbial Contaminant의 제거용
- Connection : Max id9.5 mm Tubing
- Tarsons



T521080

T521070

Open-faced Filter Holder

Cat. No.	Model	재질	Application	EA Price ₩
F03-46-113	T521070	PSF	47 mm	78,300

- Open Space 상태에서 Air Sampling • Vacuum Filtration으로 사용 가능
- Connection : Max id9.5 mm Tubing • Tarsons