


PTFE Filter Adapter for Syringe 테프론 주사기용 필터 홀더

Cat. No.	Model	사용 필터	Size	EA Price ₩
F14-460-011	N1666-08	Ø 13mm	od 21 x h35mm	435,900
F14-460-015	N1666-16	Ø 25mm	od 34 x h40mm	463,300

- 재질 : PTFE , -200 °C to 250 °C 내열
- 무게 : 14g(N1666-08), 44g(N1666-16)
- BOLA



Flow Filter / In-line Filter Holder 테프론 인라인 필터 홀더

Cat. No.	Model	사용 필터	사용 Tubing	Connection	EA Price ₩
F03-460-051	N1670-08	Ø25mm	od 3.2 /6.0mm	GL14	638,700
F03-460-055	N1670-16	Ø47mm	od 6.0/8.0mm	GL18	845,900
F03-460-059	N1670-24	Ø90mm	od 8.0/10.0mm	GL25	1,484,800


- 재질 : PTFE, PPS, -20°C to 160°C
- 용도 : Over pressure 또는 Vacuum 용
- BOLA



PFA In-line Filter / Flow Filter 테프론 인라인 필터

Cat. No.	Model	Type	사용 필터	Connecting Port	EA Price ₩
F03-460-101	N1682-08	Three Stage	Ø47mm × 3개	od 1/4"(6.35mm) Tubing	1,777,500
F03-460-105	N1678-08	Single Stage	Ø47mm × 1개		896,700


- 재질 : PFA, -200°C to 250°C
- ETFE. Glass-fiber Reinforced 연결 Nut
- Vacuum OverPressure : up to Max. 150 kPa
- 사용가능 온도 : Max 160 °C
- Three Stage : 3개의 filtering Membrane을 같이 사용할 수 있음
- BOLA



Vacuum Manifold, Stainless Steel 스테인레스 진공 메니 폴드

Cat. No.	Model	Station	EA Price ₩
F03-46-121	T50070	1-Station	521,000
F03-46-123	T50050	3-Station	1,259,000
F03-46-125	T50080	6-Station	2,431,300

- 재질 : Stainless Steel Body / PTFE Valve
- Leak Proof 2-way Valve & Vent Port
- Tarsons



Vacuum Manifold, PVC PVC 진공 메니 폴드

Cat. No.	Model	Station	EA Price ₩
F03-157-206	313400	3-Station	2,798,900
F03-157-203	313600	6-Station	4,276,300

- Non-autoclavable, up to 60°C
- 2-way Valve (PTFE Stopcock)
- Sterlitech

