

Filter Holder with Receiver, PSF PSF 리시버 일체형 필터홀더

Cat. No.	Model	Chamber 용량	Receiver 용량	EA Price ₩
F03-46-274	T50010 ①	250 ml	250 ml	107,200
F03-46-277	T50020 ②	500 ml	500 ml	140,200
F03-46-280	T50030 ③		1 L	189,800

- 재질 : PSF, -10 ~ 165°C 내열, Autoclavable
- 사용 Membrane Size : Ø47 mm
- Vacuum & Pressure Filtration용
- Graduated 상부 Chamber와 하부 Receiver
- 하부 Receiver는 2 Side Arm : id6 ~ 8 mm Vacuum Tubing 용
- Tarsons



Reusable Bottle Top Filter, GL45 PSF 바틀탑 필터, Only

Cat. No.	Model	Receiver 용량	Neck	EA Price ₩
F03-46-292	T50060	250 ml	GL 45	111,300
F03-46-298	T50061	500 ml		132,000

- 재질 : PSF, -10 ~ 165°C 내열, Autoclavable
- GL45 Thread를 가진 Bottle과 같이 사용 가능함
- Vacuum Port : PP Tubing Adapter (id6 ~ 8 mm 용)
- Glass Bottle은 별도 구매 (116 ~ 117 Page 참조)
- Tarsons



Vacuum Filters, PTFE / PFA 테프론 진공필터 리시버 일체형

- 재질 : PTFE(상부), PFA(하부), 내화학성 매우 우수
- 사용온도 : -200°C ~ +250°C
- 사용 Membrane Size : Ø47 mm / Ø90 mm
- ID 6mm 또는 8mm Vacuum Tubing 용
- 오염 물질이 표면에 정착되지 않아 세척하기 쉬우며 재사용 용이
- Vacuum & Pressure Filtration용
- BOLA

Cat. No.	Model	Diameter	Vessel 용량	Size	EA Price ₩
F03-460-150	N1650-08	Ø47 mm	240ml	od 86 x h250	3,971,000
F03-460-155	N1650-16		500ml	od 100 x h290	3,971,000
F03-460-160	N1650-24	Ø90 mm	1000ml	od 130 x h370	3,971,000

- Applications
Flow capacity for water under vacuum of 100 kPa (1000 mbar) using a PTFE filtering membrane with a thickness of 0,2 mm:
 - For membrane-Ø 47 mm with pore size 1,00 µm: Flow 510 ml/min
 - For membrane-Ø 47 mm with pore size 0,45 µm: Flow 148 ml/min
 - For membrane-Ø 47 mm with pore size 0,20 µm: Flow 57 ml/min
 - For membrane-Ø 90 mm with pore size 1,00 µm: Flow 1638 ml/min
 - For membrane-Ø 90 mm with pore size 0,45 µm: Flow 369 ml/min
 - For membrane-Ø 90 mm with pore size 0,20 µm: Flow 121 ml/min

