

Standard Filter Holder, Large 필터 홀더

Cat. No.	Model	재질	Diameter	EA Price ₩
F03-157-185	300800	SUS 304	Ø90 mm	4,008,600
F03-157-192	301900		Ø142 mm	6,551,300
F03-157-199	302500		Ø293 mm	13,534,000
F03-157-250	300900	SUS 316	Ø90 mm	5,215,300
F03-157-254	302000		Ø142 mm	8,491,000
F03-157-258	302600		Ø293 mm	17,585,500

- 재질 : Stainless Steel 304 Body
- PTFE O-ring and Silicon Rubber Gasket
- Sterlitech



Filter Holder with Reservoir 고압용 필터 홀더, Reservoir 일체형

Cat. No.	Model	Diameter	용량	재질	EA Price ₩
F03-157-822	301500	Ø47 mm	200 ml	SUS 304	5,007,700
F03-157-826	361600	Ø90 mm	750 ml		SUS 316
F03-157-830	361700			SUS 316	12,417,100
F03-157-834	302300	Ø142 mm	1500 ml	SUS 304	9,741,000
F03-157-838	302400			SUS 316	13,637,300

- Max Pressure : 5 kg/cm²
- Membrane Sealed Internally by Hand-tightening Locking Wing Nuts
- Sterlitech



Stirred Cell / High Pressure Filter Holder 고압용 필터 홀더 / 교반식 필터 홀더

Cat. No.	Model	Description	EA Price ₩
F03-157-500	HP4750	High Pressure Filter Holder	7,281,100

- 재질 : Stainless Steel 316, Buna-N Gasket and O-ring
- Connection : 1/4" FNTF (Pressure Inlet), 0.125 dia. 316 SS Tubing
- Size : w146×h199 mm (Cell wide = 51 mm)
- Volume (300 ml), Max. Pressure (1000 psig (69 bar))
- Ø47 ~ 50 mm Membrane Filter용
- High Pressure and Low Hold-up Volume (1 ml)
- Chemically Resistant Stir Cell
- Reverse Osmosis, Nanofiltration, Ultrafiltration용
- Any Standard Stirring Plate와 같이 사용 가능
- RO, NF, UF을 위한 Membrane은 425 Page 참조
- Sterlitech



구성 Component

Cat. No.	Model	Description	EA Price ₩
F03-157-610	1149782	Cell Body	4,579,200
F03-157-620	1143073	Cell Bottom	705,900
F03-157-630	1114910	Porous Support disc	258,700
F03-157-650	1143891	Cell Top	343,400
F03-157-660	1149786	High Pressure Lower Coupling	486,500
F03-157-670	1149787	High Pressure Upper Coupling	451,600
F03-157-680	1144034	Buna-N, O-ring / Gasket Kit	49,200

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
F03-157-710	1144026	Pressure Regulator Assembly, 1000 psig	2,603,900
F03-157-720	1155893	Pressure Relief Valve, 1000 psig	717,300
F03-157-730	1151427	High Pressure Hose, 6 ft	258,800
F03-157-740	1144028	PTFE O-ring / Gasket kit	98,300
F03-157-750	1144029	Viton, O-ring / Gasket kit	78,700
F03-157-760	1144027	PTFE Encapsulated Viton O-ring / Gasket kit	98,300