

Glass Vacuum Filter Holder, 25 mm 25 mm 진공 여과 장치, LukeGL®

Cat. No.	Model	Description	EA Price ₩
F03-503-100	8000015	Ø25 mm Sintered Glass Base Set	159,400
F03-503-150	8000022	Ø25 mm SUS Support Base Set	165,100

- Borosilicate Glass 3.3
- Ø25 mm Membrane Filter용
- Porosity : 40 ~ 60 um (Sintered Glass), 150 um (Stainless Steel Mesh)
- 15 ml Glass Funnel 기본 포함 • 여과병 (Filtering Flask)은 별도 구매
- LukeGL®

구성 Component, 25 mm Sintered Glass Base Set

Cat. No.	Model	Description	EA Price ₩
F03-503-510	8021521	15 ml Glass Funnel ①	45,900
F03-503-520	8021524	Aluminum Clamp, Ø25 mm ②	47,600
F03-503-530	8021522	Sintered Glass Base, Ø25 mm ③	82,900
F03-503-570	8021525	Silicon Stopper (상 Ø28 / 하 Ø25 mm)	6,900

구성 Component, 25 mm SUS Support Base Set

Cat. No.	Model	Description	EA Price ₩
F03-503-510	8021521	15 ml Glass Funnel ①	45,900
F03-503-520	8021524	Aluminum Clamp, Ø25 mm ②	47,600
F03-503-540	8021526	SUS 100 Mesh Screen, Ø25 mm ④	55,400
F03-503-550	8021527	Teflon Gasket, Ø25 mm ⑤	2,800
F03-503-560	8021523	SUS Support 용 Glass Base, Ø25 mm ⑥	50,800
F03-503-570	8021525	Silicon Stopper (상 Ø28 / 하 Ø25 mm)	6,900

Optional Accessory

Cat. No.	Model	Description	Neck id	EA Price ₩
F03-503-950	8072125	Filtering Flask, 125 ml	Ø23.5 mm	23,100
F03-503-955	8072250	Filtering Flask, 250 ml	Ø27 mm	23,400



F03-503-100



Glass Vacuum Filter Holder, 90 mm 90 mm 대용량 조인트 진공 여과 장치, LukeGL®

Cat. No.	Model	Description	EA Price ₩
F03-503-180	8000908	Ø90 mm SUS Support Base Set	680,000
F03-503-801	8000061	Ø90 mm Sintered Glass Base Set	680,000

- Borosilicate Glass 3.3
- Ø90 mm Membrane Filter용 • NS 40/35 부착형
- PTFE Coated Stainless Steel Mash Screen Set (100 mesh)
- 1 L Glass Funnel 기본 포함 • 여과병 (Filtering Flask)은 별도 구매
- LukeGL®

구성 Component, 90mm SUS Support Base Set

Cat. No.	Model	Description	EA Price ₩
F03-503-810	8934591	1 L Glass Funnel ①	210,000
F03-503-820	8934592	Aluminum Clamp, Ø90 mm	107,000
F03-503-840	8934593	SUS 100 mesh Screen, Ø90 mm ③	174,500
F03-503-850	8934595	Teflon Gasket, Ø90 mm ④	9,800
F03-503-860	8934594	NS 40/35 Glass Base ②	300,000

구성 Component, 90mm Sintered Glass Base Set

Cat. No.	Model	Description	EA Price ₩
F03-503-810	8934591	1 L Glass Funnel ①	210,000
F03-503-820	8934592	Aluminum Clamp, Ø90 mm	107,000
F03-503-804	8000063	NS 40/35 Sintered Glass Base ⑤	421,400

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
F03-503-736	8092310	Filtering Flask, 1 L, NS 40/35	65,800
F03-503-740	8092520	Filtering Flask, 2 L, NS 40/35	100,100
F03-503-960	8092740	Filtering Flask, 4 L, NS 40/35	190,800

- 안전용 플라스틱 코팅 여과병은 117 Page 참조



F03-503-801

F03-503-180

