

Glass Vacuum Filter Holder, 47 mm 47 mm 진공 여과 장치, LukeGL®

Cat. No.	Model	Description	EA Price ₩
F03-503-200	8000031	Ø47 mm Sintered Glass Base Set	205,500
F03-503-225	8000035	Ø47 mm PTFE Faced Glass Base Set	256,300
F03-503-250	8000038	Ø47 mm SUS Support Base Set	225,300

- Borosilicate Glass 3.3
- Ø47 mm Membrane Filter용
- 3 Type
 - Type 1 : Glass Fritted Support Base Set (40 ~ 60 µm)
 - Type 2 : PTFE Faced Glass Fritted Support Base Set (40 ~ 60 µm)
 - Type 3 : Stainless Steel Support Base Set (100 mesh)
- 300 ml Glass Funnel 기본 포함
- 여과병은 별도 구매
- LukeGL®

구성 Component, 47 mm Sintered Glass Base Set

Cat. No.	Model	Description	EA Price ₩
F03-503-610	8358031	300 ml Glass Funnel ①	65,100
F03-503-620	8358036	Aluminum Clamp, Ø47 mm ③	54,400
F03-503-630	8358037	Sintered Glass Base, Ø47 mm ④	97,500
F03-503-670	8358042	Silicon Stopper (상 Ø40 / 하 Ø33 mm) ⑦	10,000

구성 Component, 47 mm PTFE Faced Glass Base Set

Cat. No.	Model	Description	EA Price ₩
F03-503-613	8358032	300 ml Glass Funnel, PTFE Faced ②	128,600
F03-503-620	8358036	Aluminum Clamp, Ø47 mm ③	54,500
F03-503-633	8358044	Sintered Glass Base, PTFE Faced ⑤	117,500
F03-503-670	8358042	Silicon Stopper (상 Ø40 / 하 Ø33 mm) ⑦	10,000

구성 Component, 47 mm SUS Support Base Set

Cat. No.	Model	Description	EA Price ₩
F03-503-610	8358031	300 ml Glass Funnel ①	65,100
F03-503-620	8358036	Aluminum Clamp, Ø47 mm ③	54,500
F03-503-640	8358038	SUS 100 mesh Screen, Ø47 mm ⑧	45,900
F03-503-650	8358040	PTFE Gasket, Ø47 mm ⑨	3,000
F03-503-660	8358039	SUS Support 용 Glass Base, Ø47 mm ⑥	68,000
F03-503-670	8358042	Silicon Stopper (상 Ø40 / 하 Ø33 mm) ⑦	10,000

- 테프론이 Funnel하부와 Glass Base상부에 코팅되어져 있음
- 8358042는 1L Flask용

Optional Glass Funnel for 47 mm

Cat. No.	Model	용량	Size	EA Price ₩
F03-503-680	8358033	500 ml	Ø80×h180 mm	125,100
F03-503-690	8358034	1 L	Ø115×h180 mm	145,900

- 8358042는 1L Flask용

Filtering Flask, LukeGL® ⑩

Cat. No.	Model	용량	Size (Ø×h)	Neck id	EA Price ₩
F03-503-957	8072255	500 ml	102×185 mm	Ø31 mm	29,000
F03-503-965	8072310	1 L	125×235 mm	Ø35 mm	39,600
F03-503-970	8072320	2 L	148×275 mm	Ø38 mm	77,300

- Borosilicate Glass 3.3
- 8072320은 1420-05를 Stopper로 사용

Filtering Flask, Simax® ⑪

Cat. No.	Model	용량	Size (Ø×h)	Neck id	EA Price ₩
F04-985-096	412022250	250 ml	85×145 mm	Ø35 mm	54,100
F04-985-101	412022500	500 ml	105×175 mm		61,300
F04-985-105	412022940	1 L	135×230 mm	Ø45 mm	87,600
F04-985-109	412022950	2 L	165×255 mm	Ø60 mm	106,400

Silicon Stopper with Hole, 15 mm ⑫

Cat. No.	Model	Hole	용도	EA Price ₩
F03-132-571	1420-05	Ø15 mm	1 L, Neck id45 mm용	11,600
F03-132-573	1420-06		2 L, Neck id60 mm용	14,500

- Simax® Filtering Flask용

