

All Glass Vacuum Filter Holder, 47 mm 47 mm 죠인트 진공 여과 장치, LukeGL®

Cat. No.	Model	Description	EA Price ₩
F03-503-300	8000057	Ø47 mm Sintered Glass Base Set	301,300

- Borosilicate Glass 3.3
- Ø47 mm Membrane Filter용
- NS 40/35 Joint 부착형
- Porosity : 40 ~ 60 um (Sintered Glass)
- 300 ml Glass Funnel & 1 L Jointed Filtering Flask 기본 포함
- LukeGL®

구성 Component, 47 mm Sintered Glass Base Set

Cat. No.	Model	Description	EA Price ₩
F03-503-610	8358031	300 ml Glass Funnel ①	65,100
F03-503-620	8358036	Aluminum Clamp, Ø47 mm ②	54,500
F03-503-730	8358047	Sintered Glass Base Ø47 mm ③	173,000
F03-503-736	8092310	Filtering Flask, 1 L, NS 40/35 ④	65,800

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
F03-503-680	8358033	500 ml Glass Funnel	125,100
F03-503-690	8358034	1 L Glass Funnel	145,900
F03-503-740	8092520	Filtering Flask, 2 L, NS 40/35	100,100
F03-503-960	8092740	Filtering Flask, 4 L, NS 40/35	190,800
F03-132-750	1424-03	플라스틱 코팅 여과병, 1 L, NS 40/35	275,500
F03-132-654	1424-04	플라스틱 코팅 여과병, 2 L, NS 40/35	492,500
F03-132-658	1424-05	플라스틱 코팅 여과병, 4 L, NS 40/35	505,700



Jointed Vacuum Filter Holder, 47 mm 47 mm 죠인트 진공 여과 장치, LukeGL®

Cat. No.	Model	Description	EA Price ₩
F03-503-350	8000058	Ø47 mm SUS Support Base Set	333,500

- Borosilicate Glass 3.3
- Ø47 mm Membrane Filter용
- NS 40/35 Joint 부착형
- Stainless Steel Support Base Set (100 mesh)
- 300 ml Glass Funnel & 1 L Jointed Filtering Flask 기본 포함
- LukeGL®

구성 Component, 47 mm SUS Support Glass Base Set

Cat. No.	Model	Description	EA Price ₩
F03-503-610	8358031	300 ml Glass Funnel ①	65,100
F03-503-620	8358036	Aluminum Clamp Ø47 mm	54,500
F03-503-640	8358038	SUS 100 Mesh Screen, Ø47 mm ③	45,900
F03-503-650	8358040	PTFE Gasket, Ø47 mm ④	3,000
F03-503-733	8358049	NS 40/35 Glass Base ②	163,900
F03-503-736	8092310	Filtering Flask 1 L, NS 40/35	65,800

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
F03-503-680	8358033	500 ml Glass Funnel	125,100
F03-503-690	8358034	1 L Glass Funnel	145,900
F03-503-740	8092520	Filtering Flask, 2 L, NS 40/35	100,100
F03-503-960	8092740	Filtering Flask, 4 L, NS 40/35	190,800
F03-132-750	1424-03	플라스틱 코팅 여과병, 1 L, NS 40/35	275,500
F03-132-654	1424-04	플라스틱 코팅 여과병, 2 L, NS 40/35	492,500
F03-132-658	1424-05	플라스틱 코팅 여과병, 4 L, NS 40/35	505,700

