

Kimble Microfiltration Assembly, ULTRA-WARE® 김블 진공 여과장치, ULTRA-WARE®

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 자세한 사항은 홈페이지 (www.labmart.co.kr) 참조
- DWK Lifescience / Kimble®

25 mm Vacuum Filter Holder Set

Cat. No.	Model	Support	EA Price ₩
F03-137-400	953705-0000	Fritted Glass	463,000
F03-137-405	953730-0000	Stainless Steel	549,600

- Filtering Flask는 별도 구매

47 mm Vacuum Filter Holder Set

Cat. No.	Model	Support	EA Price ₩
F03-137-420	953755-0000	Fritted Glass	474,500
F03-137-425	953780-0000	PTFE Faced Glass	706,400
F03-137-430	953805-0000	Stainless Steel	581,700

- Filtering Flask는 별도 구매

47 mm Jointed Vacuum Filter Holder Set

Cat. No.	Model	Support	EA Price ₩
F03-137-435	953825-0000	Fritted Glass	1,129,300
F03-137-450	953855-1047	Stainless Steel	1,016,700

- 1 L Jointed Filtering Flask 포함

47 mm Vacuum Filter Holder Set with GL45 Style Bottle

Cat. No.	Model	Support	EA Price ₩
F03-137-410	953750-5347	Fritted Glass	1,027,100
F03-137-415	953800-5347	Stainless Steel	1,115,100

- 1 L Vacuum / Pressure Bottle 포함. 1.5 bar Operating Pressure

47 mm Jointed Vacuum Filter Holder Set with Solvent Pick-up Adapter

Cat. No.	Model	Support	EA Price ₩
F03-137-455	953820-1047	Fritted Glass	1,504,300

- 1 L Jointed Filtering Flask 및 od1/4" PTFE Tubing 포함

90 mm Vacuum Filter Holder Set

Cat. No.	Model	Support	EA Price ₩
F03-137-470	953805-0090	Stainless Steel	1,842,900
F03-137-475	953770-0090	Replaceable Fritted Glass	1,871,100

- Filtering Flask는 별도 구매

90 mm Jointed Vacuum Filter Holder Set

Cat. No.	Model	Description	EA Price ₩
F03-137-480	953825-0090	Fritted Glass	2,566,400
F03-137-505	953855-1090	Stainless Steel	3,238,600

- 1 L Filtering Flask 포함
- Funnel Capacity : 15 ml (25 mm), 300 ml (47 mm), 1 L (90 mm)
- Fritted Glass Support : Coarse, 40 ~ 60 um
- Stainless Steel Screen : 80 x 80 mesh

