

**Filtration Cap, ULTRA-WARE®**  
For Solvent Pickup Adaptor, Vaccum Filter Funnel  
**여과용 캡, ULTRA-WARE®**

Cat. No.	Model	Filter 규격	하부 Screw Thread	EA Price ₩
F03-137-222	953915-0047	Ø47 mm, 40 ~60 mm pore	GL 45	1,144,900

- Body 재질 : Glass-filled PTFE
- Vaccum Adaptor : 1/4" id Tubing 용
- 상부 Screw Clamp는 Solvent Pickup Adaptor 또는 Vaccum Filter System의 Funnel과 같이 사용함
- 40-60 um Fritted glass Filter 교체 가능함
- Fritted glass Filter는 PTFE/PP O-ring 포함
- DWK Lifescience / Kimble®



**Adapter Assembly with GL45 Thread, Kimble®**  
**진공여과장치용 GL45 어댑터**

Cat. No.	Model	Description	EA Price ₩
F03-137-750	179950-4532	GL45 Adapter Assembly	289,500

- Borosilicate Glass, USP Type I and ASTM E438, Type I, Class A
- Ø47 mm, Ø90 mm 진공 여과 장치용 어댑터
- 구성품 : PBT Cap (GL45, GL32), Adaptor Body, PTFE/Silicon Ring 2개, Hose Connector
- DWK Lifescience / Kimble®

**구성 Component**

Cat. No.	Model	Description	EA Price ₩
F03-137-595	73600-1413	Hose Barb Connector, 6 mm x 13-425	29,000
F03-137-598	410171-2916	PTFE/Silicon Sealing Ring, 29 mm	43,100
F03-137-600	410171-4226	PTFE/Silicon Sealing Ring, 42 mm	91,000
F03-137-603	410170-3220	PBT GL32 Open Top Cap	17,300
F03-137-605	410170-4534	PBT GL45 Open Top Cap	21,900
F03-137-755	179951-3213	Adapter Body Only	100,800

**Optional Accessory**

Cat. No.	Model	Description	EA Price ₩
F03-137-760	953907-0000	PBT Solvent Bottle Adapter, GL45, up to 4 L Bottle	279,500



- HPLC Reservoir Bottle은 118 Page 참조
- HPLC Reservoir Bottle Adapter는 119 Page 참조

