

## 2-Neck Heavy-wall Round Bottom Flask

### 2구 헤비월 환저 플라스크

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Heavy Wall Type으로 두께가 두꺼운 제품임
- 진공용으로 사용하기에도 알맞음
- Evaporating Flask로 사용 가능
- Chemglass

#### Angled Side Neck, 경사형 사이드 넥

Cat. No.	Model	용량	Joint	EA Price ₩
F24-132-171	CG-1520-41	15 ml	14/20	113,100
F24-132-174	CG-1520-43	25 ml		113,100
F24-132-177	CG-1520-46	50 ml		115,400
F24-132-180	CG-1520-49	100 ml		116,400
F24-132-183	CG-1520-02	100 ml	24/40	162,700
F24-132-186	CG-1520-04	250 ml		166,000
F24-132-189	CG-1520-10	500 ml		175,400
F24-132-192	CG-1520-14	1 L		215,100
F24-132-195	CG-1520-18	2 L		410,100

#### Vertical Side Neck, 수직형 사이드 넥

Cat. No.	Model	용량	Joint	EA Price ₩
F24-132-151	CG-1518-01	100 ml	24/40	166,400
F24-132-153	CG-1518-03	250 ml		169,200
F24-132-155	CG-1518-06	500 ml		174,900
F24-132-157	CG-1518-08	1 L		214,500
F24-132-159	CG-1518-10	2 L		402,700
F24-132-161	CG-1518-29	5 L	C=45/50, 24/40	627,600



경사형



수직형

## 3-Neck Heavy-wall Round Bottom Flask

### 3구 헤비월 환저 플라스크

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Heavy Wall Type으로 두께가 두꺼운 제품임
- 진공용으로 사용하기에도 알맞음
- Evaporating Flask로 사용 가능
- Chemglass

#### Angle Side Neck, 경사형 사이드 넥

Cat. No.	Model	용량	Joint	EA Price ₩
F24-132-301	CG-1524-A-01	15 ml	14/20	131,200
F24-132-304	CG-1524-A-02	25 ml		129,300
F24-132-307	CG-1524-A-05	50 ml		136,100
F24-132-310	CG-1524-A-08	100 ml	24/40	140,700
F24-132-313	CG-1524-01	100 ml		230,400
F24-132-316	CG-1524-82	125 ml		242,900
F24-132-319	CG-1524-02	200 ml		267,500
F24-132-322	CG-1524-03	250 ml		230,400
F24-132-325	CG-1524-38	300 ml		264,900
F24-132-328	CG-1524-05	500 ml		242,800
F24-132-331	CG-1524-07	1 L		273,600
F24-132-334	CG-1524-50	2 L		507,700
F24-132-337	CG-1524-14	3 L		561,700
F24-132-340	CG-1524-19	5 L		637,200

#### Parallel Side Neck, 수직형 사이드 넥

Cat. No.	Model	용량	Joint	EA Price ₩
F04-132-251	CG-1522-01	100 ml	24/40	248,400
F04-132-253	CG-1522-82	125 ml		256,900
F04-132-255	CG-1522-02	200 ml		253,100
F04-132-259	CG-1522-03	250 ml		229,200
F04-132-261	CG-1522-38	300 ml		253,400
F04-132-263	CG-1522-05	500 ml		241,600
F04-132-267	CG-1522-07	1 L		272,300
F04-132-271	CG-1522-50	2 L		505,300
F04-132-275	CG-1522-14	3 L		527,100
F04-132-279	CG-1522-19	5 L		635,300



경사형



수직형