


## Chemglass AIRFREE® APParatus

- Designed For The Convenient Manipulation Of Highly Reactive Air And Moisture Sensitive Compounds Using Inert Gas To Exclude Air And Moisture
- Complete Vacuum Line Assemblies And Components Using Greaseless Joint
- Heavy-wall Type으로 진공용으로 최고의 제품



### Storage Flask, AIRFREE®

#### 쉬링크 저장 플라스크, Schlenk Flask

Cat. No.	Model	용량	Joint	EA Price ₩
F24-132-561	AF-0522-01	50 ml	14/20	236,900
F24-132-563	AF-0522-02	100 ml		236,900
F24-132-567	AF-0522-03	250 ml		262,900
F24-132-571	AF-0522-04	500 ml	24/40	289,400
F24-132-565	AF-0522-06	100 ml		245,300
F24-132-569	AF-0522-07	250 ml		268,100
F24-132-573	AF-0522-05	500 ml		293,500
F24-132-575	AF-0522-08	1 L		316,300

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 50 ~ 250 ml : 0 ~ 4 mm Chem-Vac Valve
- 500 ml & 1 L : 0 ~ 8 mm Chem-Vac Valve
- Chemglass



### Storage Tube, AIRFREE®

#### 쉬링크 저장 튜브, Schlenk Tube

Cat. No.	Model	용량	Joint	EA Price ₩
F24-132-581	AF-0523-01	25 ml	14/20	231,300
F24-132-583	AF-0523-02	50 ml		267,700
F24-132-587	AF-0523-03	100 ml		301,300
F24-132-591	AF-0523-04	350 ml	24/40	332,300
F24-132-585	AF-0523-06	50 ml		277,600
F24-132-589	AF-0523-07	100 ml		310,200
F24-132-593	AF-0523-05	350 ml		344,200

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 25 ~ 50 ml : 0 ~ 4 mm Chem-Vac Valve
- 100 ml & 350 ml : 0 ~ 8 mm Chem-Vac Valve
- Chemglass



### Storage Flask, STRAUS

#### 쉬링크 저장 플라스크, Schlenk Flask

Cat. No.	Model	용량	Joint	EA Price ₩
F24-132-551	AF-0525-01	500 ml	24/40	281,000
F24-132-553	AF-0525-02	1 L		332,700

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 장기간의 Deoxygenated Solvent Storage 용
- Center Neck이 막혀 있어 Flask에 직접 Access 할 수 없음
- Unique Design permits Filling of The Flask Using 24/40 Inlet Off a Vacuum Line and Recovery of the Solvent via a Cannula or Syringe Through The Valve Body
- Chemglass

