

Quartz 1-Neck Round Bottom Flask

석영 1구 환저 플라스크, LukeGL®

Cat. No.	Model	용량	od	높이	EA Price ₩
F04-60-010	QRBF0100	100 ml	Ø65 mm	h120 mm	97,600
F04-60-025	QRBF0250	250 ml	Ø85 mm	h140 mm	156,000
F04-60-050	QRBF0500	500 ml	Ø108 mm	h170 mm	234,000
F04-60-100	QRBF1000	1 L	Ø135 mm	h215 mm	390,000
F04-60-200	QRBF2000	2 L	Ø170 mm	h250 mm	722,300

- 재질 : Quartz, 내열 온도 : 1,250°C, 상용 온도 : 1,100°C
- Joint Size : ST 24/40
- 규격외 사이즈 요청에 의해 공급 가능
- LukeGL®



PFA 1-Neck Round Bottom Flask

테프론 1구 환저 플라스크, ST 29/32

Cat. No.	Model	용량	외경 Size	EA Price ₩
F04-460-131	A158-06	100 ml	Ø 67 x 117mm	274,900
F04-460-135	A158-08	250 ml	Ø 88 x 149mm	319,500
F04-460-139	A158-09	500 ml	Ø 107 x 177mm	433,900

- 재질 : PTFE, -200 °C to 250°C 내열
- Transparent, non-porous
- Joint neck size 29
- For rotary evaporators
- BOLA



Pear-shaped Flask

Pear형 플라스크, LukeGL®

Cat. No.	Model	용량	Diameter	Socket	EA Price ₩
F04-10-288	L3031413	10 ml	Ø30 mm	14/23	9,600
F04-10-292	L3031422	25 ml	Ø40 mm		10,200
F04-10-296	L3031428	50 ml	Ø50 mm		14,300
F04-10-300	L3032937	100 ml	Ø64 mm	24/40	15,700
F04-10-304	L3032949	250 ml	Ø85 mm		21,700
F04-10-308	L3032955	500 ml	Ø105 mm		25,200

- Borosilicate Glass 3.3
- Joint Size는 요청에 의해 변경 가능함
- LukeGL®



2-Neck Pear-Shaped Flask

2구 Pear형 플라스크, LukeGL®

Cat. No.	Model	용량	Center	Side	EA Price ₩
F04-10-332	L155804	10 ml	14/23	14/23	25,300
F04-10-336	L155810	25 ml			24,900
F04-10-340	L155814	50 ml			26,000
F04-10-344	L155818	100 ml			29,900
F04-10-348	L155822	250 ml	24/40	24/40	41,300

- Borosilicate Glass 3.3
- Joint Size는 요청에 의해 변경 가능함
- LukeGL®

