

3-Neck Round Bottom Flask 3구 환저 플라스크, LukeGL®

- Borosilicate Glass 3.3
- Joint Size는 요청에 의해 변경 가능함

- LukeGL®

Angle Side-Neck, 경사형 사이드 넥

| Cat. No. | Model | 용량 | Center | Side | EA Price ₩ |
|------------|----------|--------|--------|-------|------------|
| F04-10-512 | L310104E | 50 ml | 14/23 | 14/23 | 34,300 |
| F04-10-516 | L310106E | 100 ml | | | 35,900 |
| F04-10-520 | L310106F | 100 ml | 24/40 | 24/40 | 37,200 |
| F04-10-524 | L310108F | 250 ml | | | 38,400 |
| F04-10-528 | L310110F | 500 ml | | | 45,100 |
| F04-10-532 | L310112F | 1 L | | | 50,800 |
| F04-10-536 | L310114F | 2 L | 34/45 | 24/40 | 91,900 |
| F04-10-540 | L310116F | 3 L | | | 135,500 |
| F04-10-542 | L310118F | 4 L | | | 225,600 |
| F04-10-544 | L310120F | 5 L | | | 230,500 |
| F04-10-546 | L310122F | 6 L | 45/50 | 24/40 | 322,300 |
| F04-10-548 | L310124F | 10 L | | | 504,600 |
| F04-10-552 | L310126F | 20 L | | | 695,800 |



Parallel Side-Neck, 수직형 사이드 넥

| Cat. No. | Model | 용량 | Center | Side | EA Price ₩ |
|------------|----------|--------|--------|-------|------------|
| F04-10-464 | L310204E | 50 ml | 14/23 | 14/23 | 34,700 |
| F04-10-468 | L310206E | 100 ml | | | 34,700 |
| F04-10-472 | L310206F | 100 ml | 24/40 | 24/40 | 37,200 |
| F04-10-476 | L310208F | 250 ml | | | 38,400 |
| F04-10-480 | L310210F | 500 ml | | | 45,100 |
| F04-10-484 | L310212F | 1 L | | | 50,800 |
| F04-10-488 | L310214F | 2 L | 34/45 | 24/40 | 91,900 |
| F04-10-492 | L310216F | 3 L | | | 135,500 |
| F04-10-494 | L310218F | 4 L | | | 225,600 |
| F04-10-496 | L310220F | 5 L | | | 230,500 |
| F04-10-498 | L310222F | 6 L | 45/50 | 24/40 | 322,300 |
| F04-10-500 | L310224F | 10 L | | | 504,600 |
| F04-10-504 | L310226F | 20 L | | | 695,800 |



PFA 3-Neck Round Bottom Flask 테프론 3구 환저 플라스크, ST 29/32

| Cat. No. | Model | 용량 | Side Joint | 외경 Size | EA Price ₩ |
|-------------|---------|--------|------------|---------------|------------|
| F04-460-051 | A156-12 | 100 ml | 2 x 14/23 | Ø 67 x 117mm | 1,212,900 |
| F04-460-055 | A156-20 | 250 ml | 2 x 29/32 | Ø 88 x 149mm | 2,211,800 |
| F04-460-059 | A156-36 | 500 ml | | Ø 107 x 177mm | 2,340,700 |
| F04-460-067 | A149-20 | 250 ml | 2 x GL18 | Ø 88 x 149mm | 1,199,600 |
| F04-460-071 | A149-36 | 500 ml | | Ø 107 x 177mm | 1,433,600 |

- 재질 : PTFE, -200 °C to 250°C 내열
- Center Joint : ST 29/32
- BOLA



- Chemglass, Heavy Wall Flask는 445 ~ 447 Page 참조



- Chemglass, AIRFREE® APParatus는 448 ~ 450 Page 참조

