

4-Neck Round Bottom Flask 4구 환저 플라스크, LukeGL®

- Borosilicate Glass 3.3
- Joint Size는 요청에 의해 변경 가능함
- LukeGL®

Angle Side-Neck, 경사형 사이드 넥

Cat. No.	Model	용량	Center	Side	EA Price ₩
F04-10-592	L410206E	50 ml	14/23	14/23	48,200
F04-10-596	L410208E	100 ml			48,200
F04-10-598	L410208F	250 ml	24/40	24/40	52,200
F04-10-599	L410210F	500 ml			58,800
F04-10-600	L410212F	1 L			64,600
F04-10-604	L410214F	2 L	29/42	24/40	108,200
F04-10-606	L410215F				109,100
F04-10-608	L410216F	3 L	34/45	24/40	157,300
F04-10-610	L410218F	4 L			232,200
F04-10-612	L410220F	5 L			257,600
F04-10-614	L410222F	6 L			349,400
F04-10-616	L410224F	10 L	45/50	24/40	542,200
F04-10-620	L410226F	20 L			767,500

Parallel Side-Neck, 수직형 사이드 넥

Cat. No.	Model	용량	Center	Side	EA Price ₩
F04-10-560	L410108F	250 ml	24/40	24/40	52,200
F04-10-564	L410110F	500 ml			58,800
F04-10-568	L410112F	1 L			64,600
F04-10-572	L410114F	2 L	34/45	24/40	103,700
F04-10-576	L410116F	3 L			149,800
F04-10-578	L410118F	4 L			232,200
F04-10-580	L410120F	5 L			257,600
F04-10-582	L410122F	6 L	45/50	24/40	349,400
F04-10-584	L410124F	10 L			542,200
F04-10-588	L410126F	20 L			734,900



경사형



수직형

5-Neck Round-Bottom Flask 5구 환저 플라스크, LukeGL®

- Borosilicate Glass 3.3
- Joint Size는 요청에 의해 변경 가능함
- LukeGL®

Angle Side-Neck, 경사형 사이드 넥

Cat. No.	Model	용량	Center	Side	EA Price ₩
F04-10-652	L510210F	500 ml	24/40	24/40	78,800
F04-10-656	L510212F	1 L			84,100
F04-10-660	L510214F	2 L	34/45	24/40	132,200
F04-10-664	L510216F	3 L			166,700
F04-10-666	L510218F	4 L			253,300
F04-10-668	L510220F	5 L	45/50	24/40	284,700
F04-10-670	L510222F	6 L			376,500
F04-10-672	L510224F	10 L			579,900
F04-10-676	L510226F	20 L			774,100

Parallel Side-Neck, 수직형 사이드 넥

Cat. No.	Model	용량	Center	Side	EA Price ₩
F04-10-624	L510110F	500 ml	24/40	24/40	78,800
F04-10-628	L510112F	1 L			84,100
F04-10-632	L510114F	2 L	34/45	24/40	132,200
F04-10-636	L510116F	3 L			166,700
F04-10-638	L510118F	4 L			253,300
F04-10-640	L510120F	5 L	45/50	24/40	284,700
F04-10-642	L510122F	6 L			376,500
F04-10-644	L510124F	10 L			579,900
F04-10-648	L510126F	20 L			774,100



경사형



수직형