

## Kjeldahl Flask with Plain Neck, Simax®

### 킬달 플라스크

Cat. No.	Model	용량	Size	Neck Dia.	EA Price ₩
F04-985-701	411210050	50 ml	Ø49×h200 mm	Ø22 mm	17,200
F04-985-703	411210100	100 ml	Ø60×h200 mm	Ø34 mm	22,800
F04-985-705	411210250	250 ml	Ø81×h270 mm		21,700
F04-985-707	411210500	500 ml	Ø101×h300 mm		29,800
F04-985-709	411210940	1 L	Ø126×h350 mm		35,200

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- KAVALIERGLASS

## Kjeldahl Flask with Joint

### 킬달 플라스크, LukeGL®

Cat. No.	Model	용량	Joint	Neck Dia.	EA Price ₩
F04-10-268	LFK100/3L	100 ml	24/40	Ø34 mm	37,000
F04-10-272	LFK300/3L	250 ml			46,500
F04-10-276	LFK500/3L	500 ml			42,000
F04-10-280	LFK1L/3L	1 L			44,500

- Borosilicate Glass 3.3
- 킬달 장치, 증류, 추출용
- LukeGL®

킬달 플라스크, Simax®



킬달 플라스크, LukeGL®



## Distilling Flask, Simax®

### 증류 플라스크

Cat. No.	Model	용량	Size	Neck	EA Price ₩
F04-985-751	425223100	100 ml	Ø64×h215 mm	id20 mm	33,500
F04-985-753	425223126	125 ml	Ø68×h215 mm		43,400
F04-985-755	425223250	250 ml	Ø85×h215 mm		47,300

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 고품질 내열 특급 경질 강화 유리
- 75° 각도의 Side Arm, L100 mm
- 125 ml & 250 ml는 Petroleum 계열의 증류에 사용됨
- KAVALIERGLASS



## Dissolution Flask

### 용출기 수조

Cat. No.	Model	용량	Type	Size	EA Price ₩
F04-555-201	DF1000-C173	1 L	Clear	Ø131×h173 mm	58,600

- Borosilicate Glass 3.3
- Complies with USP 711 Standard
- 50 ml 간격의 눈금
- LK Labkorea

