


## Erlenmeyer Flask, Narrow-neck, Simax®

## 삼각 플라스크

Cat. No.	Model	용량	Size (Ø×h, mm)	Neck	EA Price ₩	PK Price ₩
F04-985-025	411119025	25 ml	42×75	Ø22	15,200	10ea 142,000
F04-985-029	417119050	50 ml	51×90		15,300	10ea 143,000
F04-985-033	417119100	100 ml	64×105		17,900	10ea 167,100
F04-985-037	417119200	200 ml	79×135	Ø34	4,800	10ea 44,900
F04-985-041	417119250	250 ml	85×145		4,600	10ea 42,700
F04-985-045	417119300	300 ml	87×160		6,100	10ea 57,100
F04-985-049	417119500	500 ml	105×180	Ø42	6,800	10ea 63,100
F04-985-053	417119940	1 L	131×220		9,100	10ea 84,500
F04-985-057	417119950	2 L	166×280		11,800	6ea 65,500
F04-985-061	411119952	3 L	187×310	Ø50	52,300	6ea 292,800
F04-985-065	411119956	5 L	220×365		84,100	2ea 157,200
F04-985-069	411119968	20 L	610×370		Ø380	245,100

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Complies with ISO1773
- KAVALIERSGLASS



## Narrow Mouth Erlenmeyer Flask, Kimble®

## 세구 삼각 플라스크, White Scale

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩	PK Price ₩
F04-137-403	26500-25	25 ml	39×65 mm	8,700	12ea 96,700
F04-137-406	26500-50	50 ml	50×78 mm	10,600	12ea 118,200
F04-137-409	26500-125	125 ml	66×108 mm	8,300	12ea 92,300
F04-137-412	26500-250	250 ml	82×130 mm	10,400	12ea 116,100
F04-137-415	26500-300	300 ml	80×140 mm	11,200	12ea 125,100
F04-137-418	26500-500	500 ml	102×174 mm	10,200	6ea 57,100
F04-137-421	26500-1000	1 L	128×213 mm	19,700	6ea 110,400
F04-137-424	26500-2000	2 L	161×275 mm	42,900	4ea 160,100
F04-137-427	26500-4000	4 L	208×355 mm	101,400	-
F04-137-430	26500-6000	6 L	241×395 mm	146,400	-

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 상단 부분은 라운드 처리되어 최상의 기계적 강도를 가짐
- Body 부분은 부딪히거나 맞닿았을 때 치핑을 최소화 하기 위해 두껍게 제작됨
- 잘 지워지지 않는 Ceramic Enamel Scale
- Design Meets ASTM Specification E1404, Type I, Class I
- DWK Lifescience / Kimble®



## Narrow Mouth Erlenmeyer Flask, Kimble®

## 세구형 삼각 플라스크, Blue Scale

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩	PK Price ₩
F04-137-565	26500B-50	50 ml	50×78 mm	8,600	12ea 96,500
F04-137-568	26500B-125	125 ml	66×108 mm	7,000	12ea 78,600
F04-137-571	26500B-250	250 ml	82×130 mm	7,700	12ea 85,800
F04-137-574	26500B-500	500 ml	102×174 mm	9,700	6ea 54,300
F04-137-577	26500B-1000	1 L	128×213 mm	16,500	6ea 92,200

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 일반 플라스크보다 색상이 밝음
- Ceramic Enamel 눈금
- 눈금 색으로 구분이 가능하여 내용물 혼란의 위험을 줄여주며 식별이 용이함
- Desinged from ASTM Specification E1404, Type I, Class I Requirement
- DWK Lifescience / Kimble®



## Heavy Wall Erlenmeyer Flask, Kimble®

## 헤비월 삼각 플라스크, White Scale

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩
F04-137-448	27050-500	500 ml	104×190 mm	45,500
F04-137-451	27050-1000	1 L	135×238 mm	51,900
F04-137-454	27050-2000	2 L	168×300 mm	85,200
F04-137-457	27050-4000	4 L	208×380 mm	129,800

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 최상의 내구성을 가진 Heavy Wall Type의 플라스크
- 잘 지워지지 않는 Ceramic Enamel 눈금
- DWK Lifescience / Kimble®

