


Erlenmeyer Flask, Wide-neck, Simax®

광구 삼각 플라스크

Cat. No.	Model	용량	Size (Ø×h, mm)	Neck	EA Price ₩	PK Price ₩
F04-985-153	417106100	100 ml	64×105	Ø34	3,600	10ea 33,300
F04-985-157	417106200	200 ml	79×131	Ø50	4,800	10ea 44,900
F04-985-161	417106250	250 ml	85×140		4,100	10ea 38,200
F04-985-165	417106300	300 ml	87×156		6,100	10ea 57,100
F04-985-169	417106500	500 ml	105×175	Ø72	5,100	10ea 47,700
F04-985-173	417106940	1 L	131×220		9,300	10ea 86,500
F04-985-177	417106950	2 L	153×280		11,900	6ea 66,100

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Complies with ISO24450
- KAVALIERGLASS



Wide Mouth Erlenmeyer Flask, Kimble®

광구 삼각 플라스크, Heavy-duty

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩	PK Price ₩
F04-137-433	26650-125	125 ml	66×108 mm	6,200	12ea 68,600
F04-137-436	26650-250	250 ml	77×130 mm	9,200	12ea 102,400
F04-137-439	26650-500	500 ml	97×172 mm	14,200	6ea 79,600
F04-137-442	26650-1000	1 L	122×216 mm	24,000	6ea 134,200
F04-137-445	26650-2000	2 L	161×265 mm	46,500	4ea 173,500

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Design Meets ASTM Specification E1404, Type I, Class II
- DWK Lifescience / Kimble®



Culture Flask, Kimble®

컬처 플라스크

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩	PK Price ₩
F04-137-484	25615-125	125 ml	66×125 mm	22,700	6ea 127,000
F04-137-487	25615-250	250 ml	82×142 mm	25,300	6ea 141,700
F04-137-490	25615-500	500 ml	102×183 mm	14,500	6ea 81,000
F04-137-493	25615-1000	1 L	128×230 mm	33,100	6ea 185,300
F04-137-496	25615-2000	2 L	161×288 mm	71,000	6ea 398,100

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Neck 가열연마 처리
- Long Neck형으로 Cell Culture용 Metal Cap, Plastic Closure, Foam Plug, Cotton Plug와 함께 사용하기에 알맞음
- 잘 지워지지 않는 White Ceramic Enamel 눈금
- Neck od : Ø38 mm, 단 125 ml는 Ø25 mm임
- Design Meet ASTM Specification E1404, Type VI, Class I
- DWK Lifescience / Kimble®



Erlenmeyer Flask, PP

PP 삼각 플라스크

Cat. No.	Model	용량	Neck Size	EA Price ₩
F04-339-103	26441-0050	50 ml	19/26	2,530
F04-339-106	26441-0125	125 ml	24/29	4,100
F04-339-109	26441-0250	250 ml	29/32	6,400
F04-339-112	26441-0500	500 ml	34/35	8,300
F04-339-115	26441-1L	1 L	34/35	10,430
F04-339-118	26441-2L	2 L	29/32	14,840

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- 반투명 PP 재질, Narrow-Neck, Graduation
- LabPlasti®

