

GL45 Cell Culture Flask, Kimble®

GL45 셀 컬처 플라스크

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Blue PP GL45 Cap과 Clear Seal Ring 포함
- 입구가 넓어 내용물 주입 용이
- 잘 지워지지 않는 White Ceramic Enamel 눈금과 Marking Spot
- DWK Lifescience / Kimble®

기본형

Cat. No.	Model	용량	Cap	EA Price ₩	PK Price ₩
F04-137-505	26720-250	250 ml	GL45	80,900	6ea 453,500
F04-137-508	26720-500	500 ml		93,200	6ea 522,500
F04-137-511	26720-1000	1 L		75,200	6ea 421,600
F04-137-514	26720-2000	2 L		82,700	2ea 154,500

Safety Plastic Coated Type

Cat. No.	Model	용량	Cap	EA Price ₩	PK Price ₩
F04-137-517	KC26720-250	250 ml	GL45	104,700	6ea 587,100
F04-137-520	KC26720-500	500 ml		121,600	6ea 681,800
F04-137-523	KC26720-1000	1 L		220,500	4ea 824,300
F04-137-526	KC26720-2000	2 L		227,900	2ea 425,900

- 안전 코팅이 되어있어 유리가 깨지거나 유독성 또는 부식성 화학물질이 누출되었을 때 사용자의 위험을 줄여줌



각종 Spare Cap은 119 ~ 120 Page 참조



GL45 Cell Culture Flask, Simax®

GL45 셀 컬처 플라스크

Cat. No.	Model	용량	Size (Ø × h)	EA Price ₩	PK Price ₩
F04-985-661	8023GL/250	250 ml	85 × 150 mm	15,700	10ea 146,700
F04-985-662	8023GL/500	500 ml	105 × 105 mm	18,400	10ea 171,600
F04-985-664	8023GL/1L	1 L	135 × 225 mm	20,800	10ea 193,600
F04-985-665	8023GL/2L	2 L	160 × 290 mm	29,900	10ea 279,300

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- GL45 Blue PP Screw Cap 포함
- Heavy Wall Type의 고품질 Cell Culture용 삼각 플라스크
- KAVALIERGLASS



GL45 Baffled Cell Culture Flask, Simax®

GL45 셀 컬처 베플 플라스크

Cat. No.	Model	용량	Size (Ø × h)	EA Price ₩	PK Price ₩
F04-985-641	8023GLS/250	250 ml	85 × 150 mm	16,800	10ea 156,600
F04-985-642	8023GLS/500	500 ml	105 × 105 mm	19,500	10ea 181,600
F04-985-643	8023GLS/1L	1 L	135 × 225 mm	21,800	10ea 203,500
F04-985-644	8023GLS/2L	2 L	160 × 290 mm	29,900	10ea 279,300

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- GL45 Blue PP Screw Cap 포함
- Heavy Wall Type의 고품질 Cell Culture용 삼각 플라스크
- KAVALIERGLASS



PP Erlenmeyer Flask with Screwcap

PP 캡 삼각 플라스크

Cat. No.	Model	용량	눈금	Cap	EA Price ₩
F04-20-812	93170	50 ml	10 ml	GL45	8,100
F04-20-815	93172	100 ml	20 ml		9,400
F04-20-818	93175	250 ml	50 ml		12,700
F04-20-821	93180	500 ml	100 ml		14,700
F04-20-824	93185	1 L	200 ml		19,300

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- With Screw Cap, Graduation
- Standard Joint Stopper와 같이 사용 가능함
- Brand GmBH

