

PETG Erlenmeyer Flask, Sterile PETG 멸균 삼각 플라스크 / 컬처 플라스크

- 재질 : PETG
- Sterile : Pyrogen-free, DNase/RNase-free
- Vent Cap with 0.22 um Hydrophobic Filter : 오염 없이 기체교환 가능
- LK Labkorea

Seal Cap

Cat. No.	Model	용량	Size (mm)	Neck	CS Price ₩
F04-491-010	781001	125 ml	Ø66×h106.8	Ø34 mm	24ea 113,500
F04-491-020	782001	250 ml	Ø83×h137.8		12ea 78,700
F04-491-030	783001	500 ml	Ø101×h175.8	Ø39 mm	12ea 125,500
F04-491-040	784001	1 L	Ø127×h213.3		6ea 108,200

Vent Filter Cap

Cat. No.	Model	용량	Size (mm)	Neck	CS Price ₩
F04-491-015	781011	125 ml	Ø66×h106.8	Ø34 mm	24ea 113,500
F04-491-025	782011	250 ml	Ø83×h137.8		12ea 78,700
F04-491-035	783011	500 ml	Ø101×h175.8	Ø39 mm	12ea 125,500
F04-491-045	784011	1 L	Ø127×h213.3		6ea 108,200



PC Erlenmeyer Flask, Sterile PC 멸균 삼각 플라스크 / 컬처 플라스크

- 재질 : PC
- Sterile : Pyrogen-free, DNase/RNase-free
- Vent Cap with 0.2 um Hydrophobic Filter : 오염 없이 기체교환 가능
- 개별 멸균 비닐 포장 되어져 있어 사용이 편리함
- LK Labkorea

Flat Bottom

Cat. No.	Model	용량	Size (mm)	Neck	Cap	CS Price ₩
F04-491-050	785001	2 L	Ø162×h285	Ø47 mm	Seal	6ea 291,600
F04-491-060	786001	3 L	Ø230×h253	Ø67 mm		4ea 299,500
F04-491-055	785011	2 L	Ø162×h285	Ø47 mm	Vent Filter	6ea 291,600
F04-491-065	786011	3 L	Ø230×h253	Ø67 mm		4ea 299,500

Baffled Bottom

Cat. No.	Model	용량	Size (mm)	Neck	Cap	CS Price ₩
F04-491-110	781201	125 ml	Ø66×h107	Ø34 mm	Seal	24ea 113,500
F04-491-120	782201	250 ml	Ø83×h138			12ea 80,800
F04-491-130	783201	500 ml	Ø101×h176	Ø39 mm		12ea 125,500
F04-491-140	784201	1 L	Ø127×h213			6ea 110,600
F04-491-115	781211	125 ml	Ø66×h107	Ø34 mm	Vent Filter	24ea 113,500
F04-491-125	782211	250 ml	Ø83×h138			12ea 80,800
F04-491-135	783211	500 ml	Ø101×h176	Ø39 mm		12ea 125,500
F04-491-145	784211	1 L	Ø127×h213			6ea 110,600



Spare Cap for Flat Bottom

Cat. No.	Model	용량	Cap	CS Price ₩
F04-491-070	785903	2 L	Seal	20ea 208,400
F04-491-075	785913		Vent	20ea 226,900
F04-491-190	786903	3 L	Seal	20ea 208,400
F04-491-193	786913		Vent	20ea 226,900

High Efficiency PC Erlenmeyer Flask, Sterile 고급형 PC 멸균 삼각 플라스크 / 컬처 플라스크

- 재질 : High Transparency PC, USP Class 6 Standard, 121 °C 내열
- Sterile : Pyrogen-free, DNase/RNase-free
- Vent Cap with 0.2 um Vent Membrane : Watertight (Not Airtight)
- 관찰이 용이한 Injection Molding 눈금
- 개별 멸균 비닐 포장 되어져 있어 사용이 편리함
- LK Labkorea

Seal Cap

Cat. No.	Model	용량	Size (mm)	Neck	CS Price ₩
F04-491-150	785101	2 L	Ø162×h213.5	Ø67 mm	4ea 177,200
F04-491-160	786101	3 L	Ø162×h253.5		4ea 311,900
F04-491-163	786501			Ø90 mm	4ea 311,900
F04-491-170	787001	5 L	Ø230×h285.5		4ea 530,100

- 786501 : Wide-mouth Type

Vent Filter Cap

Cat. No.	Model	용량	Size (mm)	Neck	CS Price ₩
F04-491-155	785111	2 L	Ø162×h213.5	Ø67 mm	4ea 177,200
F04-491-165	786111	3 L	Ø162×h253.5		4ea 311,900
F04-491-172	786511			Ø90 mm	4ea 311,900
F04-491-175	787011	5 L	Ø230×h285.5		4ea 530,100

- 786511 : Wide-mouth Type



Spare Cap

Cat. No.	Model	용량	Cap	CS Price ₩
F04-491-190	786903	2 & 3 L	Seal	20ea 208,400
F04-491-193	786913		Vent	20ea 226,900
F04-491-196	787903	3 L Wide-mouth & 5 L	Seal	20ea 208,400
F04-491-199	787913		Vent	20ea 226,900