

Fernbach-type Culture Flask

기본형 컬처 플라스크

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Optional Closure는 Shaking Flask란 참조
- Chemglass

without Baffle

Cat. No.	Model	용량	Top	EA Price ₩
F04-132-511	2020-10	2.8 L	Plain Top 59 mm	323,900
F04-132-516	2020-11		Delong Top 38 mm	387,800
F04-132-526	2020-30		GL45, Vented Cap	444,300

with 3 Side Baffle

Cat. No.	Model	용량	Top	EA Price ₩
F04-132-611	2022-14	2.8 L	Plain Top 59 mm	395,900
F04-132-616	2022-16		Delong Top 38 mm	421,900
F04-132-626	2022-31		GL45, Vented Cap	493,300

with 3 Bottom Baffle

Cat. No.	Model	용량	Top	EA Price ₩
F04-132-711	2024-15	2.8 L	Plain Top 59 mm	451,800
F04-132-716	2024-17		Delong Top 38 mm	504,200
F04-132-726	2024-32		GL45, Vented Cap	510,000

with 3 Bottom & 3 Side Baffle

Cat. No.	Model	용량	Top	EA Price ₩
F04-132-811	2026-20	2.8 L	Plain Top 59 mm	564,500
F04-132-816	2026-22		Delong Top 38 mm	583,300
F04-132-826	2026-33		GL45, Vented Cap	610,100



with 3 Side Baffle



with 3 Bottom Baffle



with 3 Bottom & 3 Side Baffle



- Culture Flask용 Aluminum Stopper는 1059 Page 참조



Shaking Flask, Reinforced Top

쉐이킹 플라스크

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Reinforced Top
- Baffle의 위치 : Flask 하부 가장자리
- Chemglass

with 3 Standard Baffle Type ①

Cat. No.	Model	용량	Size	EA Price ₩
F04-132-956	2040-02	125 ml	∅70×h115 mm	70,600
F04-132-961	2040-04	250 ml	∅80×h135 mm	57,500
F04-132-971	2040-05	1 L	∅130×h215 mm	130,100
F04-132-976	2040-07	2 L	∅155×h270 mm	199,900

with 3 Deep Baffle Type ②

Cat. No.	Model	용량	Size	EA Price ₩
F04-132-861	2042-02	250 ml	∅80×h135 mm	57,100
F04-132-866	2042-04	500 ml	∅100×h180 mm	88,800
F04-132-871	2042-05	1 L	∅130×h215 mm	129,900
F04-132-876	2042-07	2 L	∅155×h270 mm	194,600

