

Shaking Flask, Delong Neck 쉐이킹 플라스크

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Delong Neck (Closure 38 mm : Option)
- Baffle의 위치 : Flask 하부 가장자리
- Chemglass

with 3 Standard Baffle Type ①

Cat. No.	Model	용량	Size	EA Price ₩
F04-132-751	2044-01	50 ml	Ø50×h105 mm	95,800
F04-132-755	2044-02	125 ml	Ø70×h140 mm	79,800
F04-132-771	2044-03	250 ml	Ø80×h160 mm	80,200
F04-132-776	2044-05	500 ml	Ø100×h200 mm	94,800
F04-132-781	2044-06	1 L	Ø130×h230 mm	108,800
F04-132-786	2044-08	2 L	Ø155×h300 mm	194,600

with 3 Deep Baffle Type ②

Cat. No.	Model	용량	Size	EA Price ₩
F04-132-551	2046-01	50 ml	Ø50×h105 mm	97,300
F04-132-556	2046-02	125 ml	Ø70×h140 mm	79,800
F04-132-561	2046-03	250 ml	Ø80×h155 mm	82,900
F04-132-566	2046-05	500 ml	Ø100×h200 mm	113,600
F04-132-571	2046-06	1 L	Ø130×h230 mm	95,900
F04-132-576	2046-08	2 L	Ø155×h300 mm	155,500
F04-132-586	2046-10	4 L	Ø205×h350 mm	416,500

with 6 Baffle Type ③

Cat. No.	Model	용량	Size	EA Price ₩
F04-132-889	2048-04	250 ml	Ø85×h160 mm	101,100
F04-132-891	2048-05	500 ml	Ø100×h200 mm	129,300
F04-132-893	2048-06	1 L	Ø130×h230 mm	162,200
F04-132-895	2048-08	2 L	Ø160×h300 mm	245,200
F04-132-897	2048-09	3 L	Ø185×h340 mm	448,000
F04-132-899	2048-10	4 L	Ø205×h350 mm	493,000

- Baffle의 위치 : Flask 하부 가장자리 & 바닥에 각 3ea

Optional Closure

Cat. No.	Model	Size	재질	EA Price ₩
F04-132-666	1502-038	38 mm	Stainless Steel ④	18,900
F04-132-674	1503-10		Polypropylene ⑤	4,600
F04-132-700	1490-038		Silicon Sponge ⑥	25,300



Shaking Flask, Screw-thread 쉐이킹 플라스크

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Baffle의 위치 : Flask 하부 가장자리
- Chemglass

with 3 Standard Baffle Type ①

Cat. No.	Model	용량	Size	EA Price ₩
F04-132-002	2051-04	250 ml	Ø80×h125 mm	187,600
F04-132-005	2051-06	500 ml	Ø100×h170 mm	228,300
F04-132-008	2051-08	1 L	Ø130×h215 mm	231,700
F04-132-011	2051-10	2 L	Ø170×h270 mm	308,800

- GL-45 PP Vented Cap (0.22 um PTFE Membrane)

with 6 Baffle Type ②

Cat. No.	Model	용량	Size	EA Price ₩
F04-132-033	2052-04	250 ml	Ø80×h125 mm	186,100
F04-132-035	2052-06	500 ml	Ø100×h170 mm	231,700
F04-132-037	2052-08	1 L	Ø130×h215 mm	239,500
F04-132-039	2052-10	2 L	Ø170×h270 mm	320,800

- GL-45 PP Vented Cap (0.22 um PVDF Membrane)

- Baffle의 위치 : Flask 하부 가장자리 & 바닥에 각 3ea

Spare Vent Cap

Cat. No.	Model	Size	Membrane	EA Price ₩
F04-132-222	1484-10	GL45	0.22 um (PTFE)	11,800

