

Roux Culture Flask, Simax® Roux 컬처 플라스크, Off-set Neck

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Neck이 Joint 타입인 제품도 공급 가능함 • KAVALLIERGLASS

Central Straight Neck ①

Cat. No.	Model	용량	Size (mm)		EA Price ₩	PK Price ₩
			높이	Neck		
F04-985-567	70/75	75 ml	125	20	20,000	10ea 186,500
F04-985-570	70/250	250 ml	180	25	22,900	16ea 341,600
F04-985-573	70/400	450 ml	195	30	26,100	8ea 194,800
F04-985-576	70/1L	1 L	255	32	44,400	6ea 248,900
F04-985-579	70/2L	2 L	320	40	63,100	5ea 294,600

Side Straight Neck ②

Cat. No.	Model	용량	Size (mm)		EA Price ₩	PK Price ₩
			높이	Neck		
F04-985-670	70A/250	250 ml	180	25	20,300	16ea 302,800
F04-985-673	70A/450	450 ml	195	30	25,700	8ea 192,000
F04-985-676	70A/1L	1 L	255	32	43,600	6ea 244,500
F04-985-679	70A/2L	2 L	320	40	62,100	5ea 290,200



Haffkine Culture Flask 하프킨 컬처 플라스크

Cat. No.	Model	용량	Size (Ø × h)	EA Price ₩
F04-555-101	CF3024H	3 L	245 × 280 mm	201,200
F04-555-104	CF4024H	4 L	246 × 310 mm	207,000

- Borosilicate Glass 3.3, ASTM E438
- Polenske Apparatus 에 적합
- LK Labkorea



Cell Culture Flask 셀 컬처 플라스크

- 재질 : High-grade Crystal-clear Polystyrene
- "TC Tested" 자체 인증 Marked
- Sterile, Pyrogen-free, Endotoxin-free, Non-cytotoxic
- DNase / RNase-free, DNA-free
- 각 제품의 Lot.NO, 유효기간 표시
- 플라스크의 모든 가장자리까지 피펫이 닿을 수 있는 구조
- Stacking 테두리로 Dish를 안전하게 쌓을 수 있음
- 0.2 um의 PTFE 필터 캡을 장착
- Sarstedt

부착 세포용 Cell Culture Flask, Color Code : Red

Cat. No.	Model	Culture Area	최적 용량	최대 용량	CS Price ₩
F04-608-001	83.3910.002	25cm ²	7 ml	12.5 ml	300ea 359,000
F04-608-004	83.3911.002	75cm ²	21 ml	55 ml	100ea 179,000
F04-608-007	83.3912.002	175cm ²	50 ml	125 ml	40ea 139,000

부유 세포용 Cell Culture Flask, Color Code : Green

Cat. No.	Model	Culture Area	최적 용량	최대 용량	CS Price ₩
F04-608-021	83.3910.502	25cm ²	7 ml	12.5 ml	300ea 501,000
F04-608-024	83.3911.502	75cm ²	21 ml	55 ml	100ea 263,000
F04-608-027	83.3912.502	175cm ²	50 ml	125 ml	40ea 209,000

