

Joint Filtering Flask

조인트 여과 플라스크, LukeGL®

Cat. No.	Model	용량	Socket	EA Price ₩
C29-11-322	FF0524	500 ml	24/40	39,800
C29-11-330	FF1024	1 L		57,800
C29-11-338	FF2024	2 L		86,700

- Borosilicate Glass 3.3
- Heavy Wall Type으로 Vacuum Application용
- 두껍고 충격에 강하며 조인트가 달려있음
- Standard Joint Socket 24/40
- LukeGL®



24/40 Joint Filtering Flask with Hose Connector

조인트 여과 플라스크, LukeGL®

Cat. No.	Model	용량	Joint	EA Price ₩
F04-503-121	FF0240	250 ml	24/40	45,800
F04-503-123	FF0540	500 ml		57,800
F04-503-125	FF1040	1 L		65,300
F04-503-127	FF2040	2 L		97,800

- Borosilicate Glass 3.3
- Heavy Wall Type으로 Vacuum Application용
- 두껍고 충격에 강하며 조인트가 달려있음
- Standard Joint Socket 24/40
- Connection : od8 mm
- LukeGL®



40/35 Joint Filtering Flask

40/35 조인트 여과 플라스크, LukeGL®

- Borosilicate Glass 3.3
- Heavy Wall Type으로 Vacuum Application용
- All Glass Vacuum Filter Holder와 같이 사용하기에 알맞음
- Standard Joint Socket 40/35
- LukeGL®

Flask Shape

Cat. No.	Model	용량	Joint	EA Price ₩
F03-503-736	8092310	1 L	40/35	65,800
F03-503-740	8092520	2 L		100,100
F03-503-960	8092740	4 L		190,800

Flask Shape, Plastic Coated, Chemglass

Cat. No.	Model	용량	Joint	EA Price ₩
F03-132-750	1424-03	1 L	40/35	275,500
F03-132-654	1424-04	2 L		492,500
F03-132-658	1424-05	4 L		505,700



PP Filtering Flask

플라스틱 여과 플라스크

Cat. No.	Model	용량	Neck	EA Price ₩
F04-20-910	26490-0250	250 ml	29/32	12,200
F04-20-920	26490-0500	500 ml	34/35	15,800
F04-20-930	26490-1L	1 L		19,800
F04-20-940	26490-2L	2 L		23,130

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Heavy Wall Type으로 고진공에 안전함
- Side Arm : od9 mm
- LabPlasti®

