

ASTM Heavy Duty Wide Mouth Volumetric Flask with Individual Serial No & Certificate, Kimble® 헤비듀 광구 메스 용량 플라스크, Class A + 개별보증서

Cat. No.	Model	용량	Tol.	Stopper	CS Price ₩
F04-137-622	92820G-5	5 ml	±0.08 ml	13	6ea 537,000
F04-137-625	92820G-10	10 ml			6ea 551,600
F04-137-628	92820G-20	20 ml			6ea 592,200
F04-137-631	92820G-25	25 ml			6ea 606,700
F04-137-634	92820G-50	50 ml			6ea 624,200
F04-137-637	92820G-100	100 ml	±0.1 ml	16	6ea 653,100
F04-137-640	92820G-200	200 ml	±0.2 ml	19	6ea 958,000
F04-137-643	92820G-250	250 ml			6ea 972,500
F04-137-646	92820G-500	500 ml			6ea 1,149,600
F04-137-649	92820G-1000	1 L	±0.3 ml	22	1ea 377,300
F04-137-652	92820G-2000	2 L	±0.05 ml	27	1ea 522,700

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 개별 Serial Number가 기재되어 있으며, 개별 보증서 포함
- 광구형으로 피펫과 사용 가능
- 헤비듀 타입으로 내구성과 안정성이 보장됨
- Glass Stopper 기본 포함
- Designed from ASTM Specification E238, Class A Serialized Requirements
- DWK Lifescience / Kimble®



ASTM Kohlrausch Volumetric Flask, Kimble® ASTM 광구 메스 용량 플라스크, Class A + 개별보증서

Cat. No.	Model	용량	Tol.	Stopper	CS Price ₩
F04-137-111	28100-100	100 ml	±0.1 ml	6.5	12ea 1,554,000
F04-137-114	28100-200	200 ml		7.5	12ea 1,724,000
F04-137-117	28100-500	500 ml	±0.2 ml	9	6ea 1,028,600

- Borosilicate Glass 3.3, USP Type I, Class A
- Calibrated to Contain (TC, In)
- Supplied with a marking spot
- Designed from ASTM Specification E438, Class A Serialized Requirements
- DWK Lifescience / Kimble®



ASTM Square Volumetric Flask, Kimble® ASTM 사각형 메스 용량 플라스크, Class A

Cat. No.	Model	용량	Tol.	Stopper	EA Price ₩
F04-137-122	28046-500	500 ml	±0.2 ml	19	295,800
F04-137-125	28046-900	900 ml	±0.3 ml	22	393,000

- Borosilicate Glass 3.3, Class A
- 유리 중량과 사각형 모양으로 인해 안정성 보장
- Supplied with a marking spot for labeling
- Glass Stopper 기본 포함
- Designed from ASTM Specification E288, Class A Serialized Requirements
- DWK Lifescience / Kimble®



ASTM Cassia Volumetric Flask, Kimble® ASTM 카시아 메스 용량 플라스크, Class A

Cat. No.	Model	용량	Tol.	Stopper	CS Price ₩
F04-137-130	28066-110	110 ml	±0.08 ml	13	6ea 696,600

- Borosilicate Glass 3.3, USP Type I, Class A
- Oils of cinnamon, chenopodium, cloves, spearmint and other essential oils 측정에 적합
- Calibrated to Contain (TC, In)
- Supplied with a marking spot
- 잘 지워지지 않는 White Ceramic Enamel 눈금
- Glass Stopper 기본 포함
- Designed from ASTM Specification E438, Class A Serialized Requirements
- DWK Lifescience / Kimble®

