

Curved Forcep, Sharp 고급형 곡형 뾰족한 핀셋

Cat. No.	Model	길이	EA Price ₩
F05-102-489	3070	115 mm	6,700
F05-102-496	3071	130 mm	7,300
F05-102-503	3072	145 mm	7,700
F05-102-510	3073	160 mm	8,700
F05-102-517	3074	200 mm	9,200

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- Usbeck



Splinter Forcep 스플린터 핀셋

Cat. No.	Model	Type	길이	EA Price ₩
F05-139-281	1310-13	Straight	130 mm	4,200
F05-139-291	1312-13	Curved		4,500

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- 앞부분이 가늘게 되어져 있는 핀셋 • 조직의 이물질 제거 및 봉합사 처리용
- LK Labkorea



Ring Forcep 링 포셉

Cat. No.	Model	Ring Size	길이	EA Price ₩
F05-695-011	91-2507	3 × id1.9 mm	110 mm	70,000
F05-695-013	91-2508	6 × id4.5 mm		70,000
F05-695-015	91-2509	3 × id1.9 mm	150 mm	81,600
F05-695-017	91-2510	6 × id4.5 mm		81,600

- 재질 : Stainless Steel 401 • 둥근 형태의 물체를 잡기에 적합함
- Sansyo



Tissue Forcep 1:2 티슈 핀셋

Cat. No.	Model	길이	EA Price ₩
F05-139-101	1304-11	115 mm	4,000
F05-139-103	1304-13	130 mm	4,200
F05-139-105	1304-14	145 mm	5,400
F05-139-107	1304-18	180 mm	8,400
F05-139-109	1304-20	200 mm	10,500
F05-139-111	1304-25	250 mm	12,900

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- 팁 끝에 예리한 고리가 있는 조직용 핀셋
- LK Labkorea



Self - Closing Forcep 셀프 클로징 핀셋

Cat. No.	Model	길이	EA Price ₩
F05-102-620	3131	105 mm	8,400

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- LK Labkorea



Soldering Forcep 납땀용 핀셋

Cat. No.	Model	Type	길이	EA Price ₩
F05-139-121	1320-12	Straight	120 mm	6,500
F05-139-123	1322-12	Curved		6,500

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- 납땀 작업시 같이 사용하기 알맞은 핀셋 • Self-closing & Reverse Action Type
- LK Labkorea

