

Tool Forcep 공구용 핀셋

Cat. No.	Model	Type	길이	EA Price ₩
F05-139-131	1324-16	Straight	160 mm	8,800
F05-139-133	1325-16	Curved		8,800

- 재질 : Stainless Steel 430
- Self-closing & Reverse Action Type
- LK Labkorea



Cover Glass Forcep 커버글라스용 핀셋, 종이 필터용 핀셋

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- LK Labkorea

Straight Type ①

길이	Stainless Steel		EA Price ₩	Teflon Coating		EA Price ₩
	Cat. No.	Model		Cat. No.	Model	
105 mm	F05-139-301	1334-10	4,000	-	-	-
115 mm	F05-102-557	3110	5,800	F05-102-571	3110T	44,000

Curved Type ②


길이	Stainless Steel		EA Price ₩	Teflon Coating		EA Price ₩
	Cat. No.	Model		Cat. No.	Model	
105 mm	F05-139-311	1336-10	4,000	-	-	-
115 mm	F05-102-564	3130	5,800	F05-102-578	3130T	44,000
130 mm	F05-102-566	15130	11,800	F05-102-576	15130T	48,700
145 mm	F05-102-568	15145	12,000	F05-102-574	15145T	49,800



Plate Forcep 플레이트용 핀셋

Cat. No.	Model	전체 길이	Tip Wide	EA Price ₩
F05-722-300	SI-300-SA	125 mm	11 mm	10,800
F05-722-310	SI-400-SA		13.5 mm	24,100
F05-722-320	SI-500-SA	130 mm	17 mm	24,100
F05-722-330	SI-600-SA		19.5 mm	24,100

- 재질 : Stainless Steel 304, 비자성 / 내부식성
- Sipel



Angular Forcep 앵글러 핀셋

Cat. No.	Model	길이	EA Price ₩
F05-139-141	1328-18	180 mm	8,700
F05-139-143	1328-20	200 mm	10,000

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- 독특한 팁의 각이진 형태의 핀셋
- LK Labkorea



Titanium Forcep 티타늄 핀셋

Cat. No.	Model	Type	길이	EA Price ₩
F05-102-663	5-TNF	Straight, Extra Sharp Tip	120mm	66,600
F05-102-667	5A-TNF	Bent, Extra Sharp Tip		66,600
F05-102-671	3C-TNF	Straight, Extra Sharp Tip	108mm	66,600
F05-102-675	7-TNF	Bent, Extra Sharp Tip		66,600

- 재질 : Titanium (6% Aluminum, 4% Vanadium), ASTM Grade 5
- 강한 내부식성, 고강도 및 Low Weight
- Blue Color Anodized
- 거의 모든 Organic Acid & Alkaline Solution과 같이 사용 가능
- Melting Point : 1,668°C
- LK Labkorea

