

Disposable IR Sample Card IR 샘플 카드

- Card Size : w5.1 × d10.2 cm
- Marking Area : 날짜, 샘플, Analyst 를 기록 할수 있음
- Universal Size 로 어떠한 Spectrophotometer 에도 사용 가능
- ICL (International Crystal Laboratory)

PTFE & Polyethylene IR Sample Card

Cat. No.	Model	Aperture	Aperture od	25ea	PK Price ₩
F11-87-008	00007362	PTFE	Ø19 mm	25ea	192,200
F11-87-016	00007365	Polyethylene		25ea	204,900

- Spectral Range : 4400 cm^{-1} to 400 cm^{-1}

Real Crystal IR Sample Card

Cat. No.	Model	Aperture	Aperture od	25ea	PK Price ₩
F11-87-035	00007095	NaCl	Ø9.5 mm	25ea	294,700
F11-87-040	00007096	KBr		25ea	363,000
F11-87-045	00007097	KCl		25ea	313,700
F11-87-050	00007131	NaCl	Ø15 mm	25ea	388,400
F11-87-055	00007132	KBr		25ea	388,400
F11-87-060	00007133	KCl		25ea	394,800
F11-87-065	00007087	NaCl	Ø19 mm	25ea	363,000
F11-87-070	00007090	KBr		25ea	429,200
F11-87-075	00007091	KCl		25ea	377,800

- Non-absorbing Substrate의 IR 분석시 이용됨
- Liquid, Mull, Film casting, Pyrolysis Sample용으로 적합
- PTFE, Polyethylene Sample Card와 달리 흡광도 Peak가 나타나지 않음

Cover Slip

Cat. No.	Model	Crystal	Crystal od	25ea	PK Price ₩
F11-87-110	00007092	NaCl	Ø19 mm	25ea	251,800
F11-87-115	00007093	KBr		25ea	294,700
F11-87-120	00007094	KCl		25ea	269,000

- Crystal Cover Slip : 시료가 가루이거나 휘발성이 있는 액체일때 뒷개로 이용,
Ø19 mm Aperture IR Card용임



Unpolished Crystal Optic 크리스탈 윈도우

Cat. No.	Model	Crystal	Crystal Size	EA Price ₩
F11-87-150	0002L589	NaCl	Ø13 × h2 mm	44,800
F11-87-154	0002L131		Ø25 × h5 mm	42,700
F11-87-158	0002L369		Ø32 × h3 mm	71,700
F11-87-162	0002L4051		Ø32 × h3 mm, drilled	81,100
F11-87-166	0002M553	KCl	Ø13 × h2 mm	58,900
F11-87-170	0002M249		Ø25 × h5 mm	44,800
F11-87-174	0002M411		Ø32 × h3 mm	74,700
F11-87-178	0002M3075		Ø32 × h3 mm, drilled	83,400
F11-87-182	0002N396	KBr	Ø13 × h2 mm	51,300
F11-87-186	0002N88		Ø25 × h5 mm	48,100
F11-87-190	0002N90		Ø32 × h3 mm	80,100
F11-87-194	0002N534		Ø32 × h3 mm, drilled	85,300

- 완전 평면 Window임
- 그 외 사이즈는 문의 바람
- ICL (International Crystal Laboratory)

