

## Glass Filter Funnel

### 깔때기형 글라스 필터, ISO

| Cat. No.    | Model | 용량     | Porosity | 필터 (Ø) | Stem (Ø) | Size (Ø × h) | EA Price ₩ |
|-------------|-------|--------|----------|--------|----------|--------------|------------|
| F07-228-004 | 21501 | 50 ml  | P1       | 36 mm  | 10 mm    | 40×53 mm     | 46,600     |
| F07-228-008 | 21502 |        | P2       |        |          |              | 46,600     |
| F07-228-012 | 21503 |        | P3       |        |          |              | 46,600     |
| F07-228-016 | 21504 |        | P4       |        |          |              | 50,100     |
| F07-228-020 | 21505 |        | P5       |        |          |              | 85,400     |
| F07-228-024 | 21751 | 75 ml  | P1       | 48 mm  | 10 mm    | 52×50 mm     | 52,400     |
| F07-228-028 | 21752 |        | P2       |        |          |              | 52,400     |
| F07-228-032 | 21753 |        | P3       |        |          |              | 52,400     |
| F07-228-036 | 21754 |        | P4       |        |          |              | 57,800     |
| F07-228-040 | 21755 |        | P5       |        |          |              | 94,200     |
| F07-228-044 | 21121 | 125 ml | P1       | 60 mm  | 10 mm    | 65×56 mm     | 73,700     |
| F07-228-048 | 21122 |        | P2       |        |          |              | 73,700     |
| F07-228-052 | 21123 |        | P3       |        |          |              | 73,700     |
| F07-228-056 | 21124 |        | P4       |        |          |              | 81,900     |
| F07-228-060 | 21125 |        | P5       |        |          |              | 146,600    |
| F07-228-064 | 21251 | 250 ml | P1       | 80 mm  | 18 mm    | 85×65 mm     | 110,500    |
| F07-228-068 | 21252 |        | P2       |        |          |              | 110,500    |
| F07-228-072 | 21253 |        | P3       |        |          |              | 110,500    |
| F07-228-076 | 21254 |        | P4       |        |          |              | 121,800    |
| F07-228-080 | 21255 |        | P5       |        |          |              | 185,400    |
| F07-228-084 | 21551 | 500 ml | P1       | 95 mm  | 22 mm    | 100×100 mm   | 135,400    |
| F07-228-088 | 21552 |        | P2       |        |          |              | 135,400    |
| F07-228-092 | 21553 |        | P3       |        |          |              | 135,400    |
| F07-228-096 | 21554 |        | P4       |        |          |              | 161,100    |
| F07-228-100 | 21555 |        | P5       |        |          |              | 314,300    |
| F07-228-104 | 21111 | 1 L    | P1       | 124 mm | 22 mm    | 130×115 mm   | 219,800    |
| F07-228-108 | 21112 |        | P2       |        |          |              | 219,800    |
| F07-228-112 | 21113 |        | P3       |        |          |              | 219,800    |
| F07-228-116 | 21114 |        | P4       |        |          |              | 287,800    |
| F07-228-120 | 21115 |        | P5       |        |          |              | 530,600    |
| F07-228-124 | 21401 | 4 L    | P1       | 180 mm | 30 mm    | 190×206 mm   | 714,300    |
| F07-228-128 | 21402 |        | P2       |        |          |              | 714,300    |
| F07-228-132 | 21403 |        | P3       |        |          |              | 714,300    |
| F07-228-136 | 21404 |        | P4       |        |          |              | 824,800    |
| F07-228-144 | 21801 | 8 L    | P1       | 205 mm | 40 mm    | 215×280 mm   | 1,408,300  |
| F07-228-148 | 21802 |        | P2       |        |          |              | 1,408,300  |
| F07-228-152 | 21803 |        | P3       |        |          |              | 1,408,300  |
| F07-228-156 | 21804 |        | P4       |        |          |              | 1,626,500  |

- Borosilicate Glass 3.3
- ISO 규격에 따른 Porosity
- Heavy Wall Type으로 VitraPOR® 고품질 Filter 사용
- Robu / VitraPOR®



### ISO 규격에 따른 Porosity

| Porosity | ISO4793 | Pore Size    |
|----------|---------|--------------|
| 1        | P160    | 100 ~ 600 um |
| 2        | P100    | 40 ~ 100 um  |
| 3        | P40     | 16 ~ 40 um   |
| 4        | P16     | 10 ~ 16 um   |
| 5        | P1.6    | 1.0 ~ 1.6 um |

|   |  |
|---|--|
| 1 | Coarse filtration, gas distribution, liquid distribution, gas cleaning   |
| 2 | Preparative fine filtration, mercury filtration, preparative work with coarse crystalline products               |
| 3 | Analytical work with medium fine precipitates analytical filtration, fine gas filtration, special filtration     |
| 4 | Analytical fine filtration, analytical and preparative works with the finest precipitates return and stop valves |
| 5 | General bacteria filtration, sterile filtration  |

## Glass Funnel with Filter

### 필터 부착형 깔때기

| Cat. No.   | Model | Porosity | 용량     | 필터 (Ø) | Stem (Ø) | EA Price ₩ |
|------------|-------|----------|--------|--------|----------|------------|
| F07-83-366 | 22261 | 1        | 170 ml | 30 mm  | 10 mm    | 69,600     |
| F07-83-372 | 22262 | 2        |        |        |          | 69,600     |
| F07-83-378 | 22263 | 3        |        |        |          | 69,600     |
| F07-83-384 | 22264 | 4        |        |        |          | 75,500     |
| F07-83-390 | 22265 | 5        |        |        |          | 106,700    |

- Borosilicate Glass 3.3
- ISO 규격에 따른 Porosity
- Heavy Wall Glass
- Robu / VitraPOR®

