

Glass Filter Funnel 깔때기형 글라스 필터, ASTM

Cat. No.	Model	용량	Porosity	필터 (Ø)	Stem (Ø)	Height	EA Price ₩
F07-132-389	1402-09	30 ml	XC	30 mm	8 mm	45 mm	132,200
F07-132-393	1402-10		C				95,900
F07-132-397	1402-11		M				95,900
F07-132-401	1402-12		F				95,900
F07-132-405	1402-13	60 ml	XC	40 mm	10 mm	50 mm	108,600
F07-132-409	1402-14		C				108,600
F07-132-413	1402-15		M				108,600
F07-132-417	1402-16		F				108,600
F07-132-421	1402-18	150 ml	C	60 mm	16 mm	60 mm	191,800
F07-132-425	1402-19		M				191,800
F07-132-429	1402-20		F				191,800
F07-132-433	1402-22		C				269,800
F07-132-437	1402-23	350 ml	M	80 mm	19 mm	75 mm	269,800
F07-132-451	1402-24		F				269,800
F07-132-455	1402-25		XC				323,900
F07-132-459	1402-26	600 ml	C	90 mm	25 mm	100 mm	323,900
F07-132-463	1402-27		M				323,900
F07-132-467	1402-28		F				323,900
F07-132-471	1402-30		C				706,000
F07-132-475	1402-31	2 L	M	125 mm	25 mm	145 mm	706,000
F07-132-479	1402-32		F				857,500
F07-132-483	1402-34		C				920,200
F07-132-487	1402-35	3 L	M	150 mm	25 mm	175 mm	920,200
F07-132-491	1402-36		F				1,014,000

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- ASTM F128-94 규격에 따른 Porosity
- Heavy Wall Type으로 사용이 용이함
- Height : 필터 Disc 위로의 높이임
- Chemglass



ASTM F128-94 규격에 따른 Porosity

Porosity Grade	ASTM Grade	Pore Size
Fine	F	4.0 ~ 5.5 um
Medium	M	10 ~ 16 um
Coarse	C	40 ~ 60 um
Extra Coarse	XC	170 ~ 220 um

F	Filtration and Extraction
M	
C	Gas Dispersion, Gas Washing and Absorption. Coarse Filtration.
XC	

Glass Filter Funnel 깔때기형 글라스 필터, LukeGL®

Cat. No.	Model	용량	Porosity	필터 (Ø)	EA Price ₩
F07-503-100	FNG030G1	35 ml	G1	30 mm	12,300
F07-503-103	FNG030G2		G2		12,300
F07-503-106	FNG030G3		G3		12,300
F07-503-109	FNG030G4		G4		12,300
F07-503-112	FNG040G1	80 ml	G1	40 mm	27,500
F07-503-115	FNG040G2		G2		27,500
F07-503-118	FNG040G3		G3		27,500
F07-503-121	FNG040G4		G4		27,500
F07-503-124	FNG065G1	200 ml	G1	65 mm	46,600
F07-503-127	FNG065G2		G2		46,600
F07-503-130	FNG065G3		G3		46,600
F07-503-133	FNG065G4		G4		46,600
F07-503-136	FNG090G1	500 ml	G1	90 mm	92,500
F07-503-139	FNG090G2		G2		92,500
F07-503-142	FNG090G3		G3		92,500
F07-503-145	FNG090G4		G4		92,500
F07-503-148	FNG120G1	1 L	G1	120 mm	251,400
F07-503-151	FNG120G2		G2		251,400
F07-503-154	FNG120G3		G3		251,400
F07-503-157	FNG120G4		G4		251,400

- Borosilicate Glass 3.3
- LukeGL®



Funnel Porosity

Porosity Grade	Grade	Pore Size
Fine	G4	5 ~ 15 um
Medium	G3	15 ~ 40 um
Coarse	G2	40 ~ 90 um
Extra Coarse	G1	90 ~ 150 um