


Glass Micro-filter Funnel 마이크로 필터 깔때기, LukeGL®

Cat. No.	Model	Porosity	용량	Disc dia.	Stem	EA Price ₩
F07-52-926	L378002M	1	8 ml	Ø15 mm	Ø8 mm	21,700
F07-52-932	L378004M	2				21,700
F07-52-938	L378006M	3				25,300

- Borosilicate Glass 3.3
- ISO 규격에 따른 Porosity
- 요청에 의해 규격외 제품도 공급 가능
- 길이 : 135 mm
- LukeGL®



Pipeline Filter 파이프라인 필터

Cat. No.	Model	Porosity	필터 (Ø)	Stem (Ø)	높이	EA Price ₩
F07-228-621	23301	P1	30 mm	12 mm	140 mm	69,900
F07-228-623	23302	P2				69,900
F07-228-625	23303	P3				69,900
F07-228-627	23304	P4				76,300
F07-228-629	23305	P5				110,000
F07-228-631	23601	P1	60 mm	14 mm	230 mm	158,400
F07-228-633	23602	P2				158,400
F07-228-635	23603	P3				158,400
F07-228-637	23604	P4				169,500
F07-228-639	23605	P5				217,500
F07-228-641	23901	P1	90 mm	18 mm	280 mm	285,700
F07-228-643	23902	P2				285,700
F07-228-645	23903	P3				285,700
F07-228-647	23904	P4				307,600
F07-228-649	23905	P5				433,400

- Borosilicate Glass 3.3
- ISO 규격에 따른 Porosity
- Heavy Wall Type으로 VitraPOR® 고품질 Filter 사용
- Robu / VitraPOR®



Filter Adapter 필터 어댑터

- LK Labkorea

Neoprene Rubber Filter Adapter Set

Cat. No.	Model	Description	EA Price ₩
F07-11-010	L140101	Neoprene Filter Adapter Set, 7종 세트 ①	8,300
F07-11-011	L140102	Large Filter Adapter ②	16,900

- L140101 : 7개의 서로 다른 크기의 어댑터를 Set로 구성 (No.1 ~ 7)
- Neoprene 재질로 Silicon보다 딱딱함

Silicon Filter Adapter Set

Cat. No.	Model	Description	EA Price ₩
F07-11-014	L140111	Silicon Filter Adapter Set, 6종 세트 ③	9,700

- 6개의 서로 다른 경사도의 어댑터를 Set로 구성 (No.1 ~ 6)



Neoprene Rubber



Neoprene Filter Adapter



Silicon

Silicon

No	Top od	Bottom od	Thickness
1	Ø21 mm	Ø11 mm	2t
2	Ø25 mm	Ø15 mm	2t
3	Ø37 mm	Ø22 mm	3t
4	Ø45 mm	Ø28 mm	3t
5	Ø55 mm	Ø36 mm	3t
6	Ø68 mm	Ø46 mm	4t

Neoprene

No	Top od	Bottom od	Thickness
1	Ø21 mm	Ø11 mm	2t
2	Ø27 mm	Ø16 mm	3t
3	Ø37 mm	Ø22 mm	3t
4	Ø47 mm	Ø29 mm	3t
5	Ø58 mm	Ø38 mm	4t
6	Ø68 mm	Ø45 mm	5t
7	Ø86 mm	Ø57 mm	5t
Large	Ø97 mm	Ø67 mm	6t