

Plastic Disposable Filter Funnel, Barrel Shape 대용량 일회용 플라스틱 필터 펀넬, 2L

Cat. No.	Model	용량	Porosity	Description	PK Price ₩
F07-132-181	OP-6602-J-20	2 L	20 um	Funnel only	8ea 637,800
F07-132-183	OP-6602-J-21			Funnel with Lid	8ea 671,900

- 재질 : PP Body, 20 Micron PE Fritted Disc, Stem Length : 70 mm
- Molded PP Funnel with Perforated Support Plate
- Vacuum Adapter(CG-1049-M, 별도 구매)와 같이 사용하기에 알맞음
- Chemglass

Optional Vacuum Adapter

Cat. No.	Model	상부 Outer Joint	하부 Inner Joint	Connection	EA Price ₩
F07-132-187	CG-1049-M-24	24/25	24/25	od10 mm	139,900
F07-132-189	CG-1049-M-45		45/50		186,300



Plastic Disposable Filter Funnel 일회용 플라스틱 필터 깔때기

- 재질 : PP Body, Stem Length : 80 mm
- 10 Micron PE Fritted Disc, 단 OP-6602-16 & 18은 20 Micron / OP-6602-13은 40 Micron
- For Solid-Liquid Separation in Chemistry & Biological
- Vacuum Adapter(CG-1049-J, 별도구매)와 같이 사용하기에 알맞음
- Chemglass

Shelf Pack in Convenience Tray

Cat. No.	Model	용량	Porosity	Funnel Top	PK Price ₩
F07-132-151	OP-6602-C10	18 ml	10 um	Flared	24ea 187,200
F07-132-153	OP-6602-C12	40 ~ 60 ml			24ea 208,100
F07-132-155	OP-6602-C14	110 ml		Straight	12ea 163,600

Bulk Pack in Vinly Zipper Bag

Cat. No.	Model	용량	Porosity	Funnel Top	PK Price ₩
F07-132-161	OP-6602-10	18 ml	10 um	Flared	100ea 646,400
F07-132-163	OP-6602-12	40 ~ 60 ml		Straight	100ea 708,800
F07-132-164	OP-6602-13		40 um	Flared	100ea 714,400
F07-132-165	OP-6602-14	110 ml	10 um	Straight	50ea 541,400
F07-132-166	OP-6602-16	300 ml	20 um		30ea 606,700
F07-132-167	OP-6602-18	500 ml			20ea 569,900



Plastic Disposable Filter Funnel, Cone Shape 일회용 플라스틱 필터 깔때기

Cat. No.	Model	용량	Porosity	Funnel Top	PK Price ₩
F07-132-170	OP-6605-08	8 ml	10 um	Cone Shape	100ea 538,700
F07-132-171	OP-6605-10	20 ml			100ea 654,800
F07-132-173	OP-6605-12	60 ml			100ea 717,200

- 재질 : PP Body, 10 Micron PE Fritted Disc, Stem Length : 80 mm
- Molded PP Funnel with Perforated Support Plate
- Vacuum Adapter(CG-1049-J, 별도구매)와 같이 사용하기에 알맞음
- Chemglass



Vacuum Adapter with Stopcock 조절식 진공 어댑터

Cat. No.	Model	Bottom Connection	Vacuum	EA Price ₩
F07-132-141	CG-1049-J-14	14/20	Bore 2 mm, od10 mm Hose Nipple	246,700
F07-132-143	CG-1049-J-24	24/25		254,200
F07-132-145	CG-1049-J-29	29/26		254,200
F07-132-147	CG-1049-J-30	20-400 GPI Thread		247,100

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 일회용 플라스틱 필터 Funnel과 사용하기에 적합
- Chemglass

