

Glass Filter Funnel with Vacuum Adapter and Inner Joint 24/40 진공 여과 유리 깔때기, ISO

Cat. No.	Model	용량	Porosity	필터 (Ø)	EA Price ₩
F07-228-416	21751/3S	75 ml	1	Ø48 mm	79,200
F07-228-420	21752/3S		2		79,200
F07-228-424	21753/3S		3		79,200
F07-228-428	21754/3S		4		82,200
F07-228-432	21755/3S		5		126,500
F07-228-436	21121/3S	125 ml	1	Ø60 mm	94,300
F07-228-440	21122/3S		2		94,300
F07-228-444	21123/3S		3		94,300
F07-228-448	21124/3S		4		107,000
F07-228-452	21125/3S		5		187,500
F07-228-456	21251/3S	250 ml	1	Ø80 mm	130,400
F07-228-460	21252/3S		2		130,400
F07-228-464	21253/3S		3		130,400
F07-228-468	21254/3S		4		158,200
F07-228-472	21255/3S		5		225,300
F07-228-476	21551/3S	500 ml	1	Ø95 mm	176,400
F07-228-480	21552/3S		2		176,400
F07-228-484	21553/3S		3		176,400
F07-228-488	21554/3S		4		195,600
F07-228-492	21555/3S		5		376,700
F07-228-496	21111/3S	1 L	1	Ø124 mm	256,500
F07-228-500	21112/3S		2		256,500
F07-228-504	21113/3S		3		256,500
F07-228-508	21114/3S		4		321,700
F07-228-512	21115/3S		5		623,100

- Borosilicate Glass 3.3
- ISO 규격에 따른 Porosity
- Heavy Wall Type으로 VitraPOR® 고품질 Filter 사용
- ST24/40, od10 mm Connection
- Robu / VitraPOR®



ISO 규격에 따른 Porosity

Porosity	ISO4793	Pore Size
0	P250	160 ~ 250 um
1	P160	100 ~ 600 um
2	P100	40 ~ 100 um
3	P40	16 ~ 40 um
4	P16	10 ~ 16 um
5	P1.6	1.0 ~ 1.6 um

Glass Filter Funnel with Vacuum Adapter and Inner Joint 24/40 진공 여과 유리 깔때기, ASTM

Cat. No.	Model	용량	Porosity	필터 (Ø)	EA Price ₩
F07-132-293	1406-04	30 ml	C	30 mm	212,000
F07-132-297	1406-05		M		212,000
F07-132-301	1406-06		F		212,000
F07-132-305	1406-07	60 ml	C	40 mm	274,900
F07-132-309	1406-08		M		274,900
F07-132-313	1406-09		F		274,900
F07-132-317	1406-10	150 ml	C	60 mm	363,400
F07-132-321	1406-11		M		363,400
F07-132-325	1406-12		F		363,400
F07-132-329	1406-24	350 ml	C	80 mm	418,800
F07-132-333	1406-25		M		418,800
F07-132-337	1406-26		F		423,700
F07-132-341	1406-30	600 ml	C	90 mm	506,500
F07-132-345	1406-31		M		506,500
F07-132-349	1406-32		F		506,500

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- ASTM 규격에 따른 Porosity
- ST24/40, od10 mm Connection
- Chemglass



ASTM 규격에 따른 Porosity

Porosity Grade	ASTM Grade	Pore Size
Fine	F	4.0 ~ 5.5 um
Medium	M	10 ~ 16 um
Coarse	C	40 ~ 60 um

F	Filtration and Extraction
M	
C	Gas Dispersion, Gas Washing and Absorption. Coarse Filtration.