

## Glass Filter Funnel with Vacuum Adapter and Inner Joint 24/40 진공 여과 깔때기, LukeGL®

Cat. No.	Model	용량	Porosity	필터 (Ø)	EA Price ₩
F07-503-201	FVG030G1	35 ml	G1	30 mm	41,000
F07-503-204	FVG030G2		G2		41,000
F07-503-207	FVG030G3		G3		41,000
F07-503-210	FVG030G4		G4		43,000
F07-503-213	FVG040G1	80 ml	G1	40 mm	56,000
F07-503-216	FVG040G2		G2		56,000
F07-503-219	FVG040G3		G3		56,000
F07-503-222	FVG040G4		G4		58,000
F07-503-225	FVG065G1	200 ml	G1	65 mm	77,000
F07-503-228	FVG065G2		G2		77,000
F07-503-231	FVG065G3		G3		77,000
F07-503-234	FVG065G4		G4		79,000
F07-503-237	FVG090G1	500 ml	G1	90 mm	126,500
F07-503-240	FVG090G2		G2		126,500
F07-503-243	FVG090G3		G3		126,500
F07-503-246	FVG090G4		G4		135,500
F07-503-249	FVG120G1	1 L	G1	120 mm	289,200
F07-503-252	FVG120G2		G2		289,200
F07-503-255	FVG120G3		G3		289,200
F07-503-258	FVG120G4		G4		291,000

- Borosilicate Glass 3.3
- LukeGL®



### Funnel Porosity

Porosity Grade	Grade	Pore Size
Fine	G4	5 ~ 15 um
Medium	G3	15 ~ 40 um
Coarse	G2	40 ~ 90 um
Extra Coarse	G1	90 ~ 150 um

## Cylindrical Glass Filter Funnel with Vacuum Adapter and Inner Joint 진공 여과 유리 깔때기, Column Type, LukeGL®

Cat. No.	Model	용량	Pore	od	Height	EA Price ₩
F07-10-144	L140661	450 ml	2	Ø70 mm	150 mm	106,000
F07-10-150	L140662		3			106,000
F07-10-156	L140663		4			108,400
F07-10-162	L140664		2			108,400
F07-10-168	L140665	600 ml	3	Ø70 mm	200 mm	108,400
F07-10-174	L140666		4			106,300
F07-10-180	L140667		2			114,500
F07-10-186	L140668		3			114,500
F07-10-192	L140669	700 ml	4	Ø70 mm	250 mm	118,100
F07-10-202	L140670		2			118,100
F07-10-208	L140671		3			118,100
F07-10-214	L140672		4			126,500
F07-10-220	L140673	850 ml	2	Ø90 mm	150 mm	150,600
F07-10-226	L140674		3			150,600
F07-10-232	L140675		4			162,700
F07-10-238	L140676		2			156,600
F07-10-244	L140677	1.2 L	3	Ø90 mm	200 mm	156,600
F07-10-250	L140678		4			157,500
F07-10-256	L140679	1.5 L	2	Ø90 mm	250 mm	157,500
F07-10-262	L140680		3			157,500
F07-10-268	L140681		4			174,300
F07-10-274	L140682		2			168,700
F07-10-280	L140683	1.7 L	3	Ø90 mm	300 mm	168,700
F07-10-286	L140684		4			185,500

- Borosilicate Glass 3.3
- 기본 Joint Size : 24/40
- Height : Filter Disk 위로의 높이
- 요청에 의해 규격의 사이즈 공급 가능
- ISO 규격에 따른 Porosity
- LukeGL®

