

PTFE Vacuum Filter Funnel with Joint 테프론 진공 여과 깔때기

Cat. No.	Model	Joining Size	사용 필터	용량	EA Price ₩
F07-460-011	N1658-08	29/32	Ø 47mm	125 ml	1,238,200

- 재질 : PFA Chamber , -200 °C to 250°C 내열
PTFE Filtration Unit , -200 °C to 250°C 내열
- Hose Connection : id 6mm~8mm
- Size : od62 x h188mm
- BOLA

Optional Reduction Adapter

Cat. No.	Model	Description	규격	EA Price ₩
A01-83-068	1502924	Reduction Adapter	29/32 to 24/40	21,400



Jacketed Filter Funnel 자켓식 진공 여과 깔때기

Cat. No.	Model	용량	Porosity	필터 (Ø)	EA Price ₩
F07-132-619	CG-1403-02	30 ml	C	30 mm	562,000
F07-132-623	CG-1403-03		M		562,000
F07-132-627	CG-1403-06	60 ml	C	40 mm	619,600
F07-132-631	CG-1403-07		M		619,600
F07-132-635	CG-1403-10	150 ml	C	60 mm	663,400
F07-132-639	CG-1403-11		M		663,400
F07-132-643	CG-1403-14	350 ml	C	80 mm	858,100
F07-132-647	CG-1403-15		M		873,700

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Stem Length : 60 mm (Jacket 밑 부분의 길이)
- Chemglass



Pressure Filter Funnel 정압식 필터 깔때기

Cat. No.	Model	용량	Disc od	EA Price ₩
F07-132-701	CG-1405-01	150 ml	50 mm	589,300
F07-132-703	CG-1405-03	250 ml	60 mm	589,300
F07-132-705	CG-1405-05	500 ml	70 mm	839,500
F07-132-707	CG-1405-07	1 L	80 mm	991,200

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Porosity : Coarse, 40 ~ 60 um
- Positive Pressure Filtration용, up to 15 psi
- O-ring Joint : Large, Greasefree Opening
- 구성 : Filter Funnel, Top Inlet Adapter, Viton O-ring, Pinch Clamp
- Chemglass

Filter Funnel Only


Cat. No.	Model	용량	O-ring Joint	Clamp	EA Price ₩
F07-132-711	CG-1405-10	150 ml	#30	#50	371,900
F07-132-713	CG-1405-12	250 ml	#40	#65	432,300
F07-132-715	CG-1405-14	500 ml	#50	#75	498,600
F07-132-717	CG-1405-16	1 L			622,400

Top Adapter Only

Cat. No.	Model	O-ring Joint Size	EA Price ₩
F07-132-721	CG-1405-20	#30	127,200
F07-132-723	CG-1405-21	#40	137,800
F07-132-725	CG-1405-22	#50	168,100

Spare Pinch Clamp

Cat. No.	Model	O-ring Joint Size	EA Price ₩
V08-132-626	CG-150-08	#30	73,600
V08-132-624	CG-150-09	#40	102,500
V08-132-620	CG-150-21	#50	136,300



Spare Viton O-ring

Cat. No.	Model	O-ring Joint Size	PK Price ₩
V08-132-912	CG-305-223	#30	3ea 18,600
V08-132-913	CG-305-226	#40	3ea 21,700
P05-132-579	CG-305-229	#50	3ea 25,900