

Glass Buchner Funnel 부흐너 유리 깔때기

Cat. No.	Model	용량	Size (Ø×h)	적용 필터	EA Price ₩
F07-53-009	213412207	70 ml	48×132 mm	Ø45 mm	112,500
F07-53-018	213412807	125 ml	60×140 mm	Ø55 mm	132,100
F07-53-027	213413409	220 ml	73×190 mm	Ø70 mm	223,000
F07-53-036	213414405	500 ml	95×240 mm	Ø90 mm	225,300
F07-53-045	213415401	1 L	120×270 mm	Ø110 mm	332,700

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Filter-disc가 아닌 많은 구멍이 뚫려 있는 Perforated-plate
- DWK Lifescience / Lenz



Glass Buchner Funnel with Vacuum Adapter and Inner Joint 24/40 진공 여과 부흐너 유리 깔때기, LukeGL®

Cat. No.	Model	용량	Size (Ø×h)	적용 필터	EA Price ₩
F07-10-009	L140691	70 ml	48×130 mm	Ø45 mm	152,800
F07-10-018	L140692	125 ml	60×140 mm	Ø55 mm	180,000
F07-10-027	L140693	220 ml	73×180 mm	Ø70 mm	255,900
F07-10-036	L140694	500 ml	95×240 mm	Ø90 mm	261,800
F07-10-045	L140695	1 L	120×280 mm	Ø110 mm	375,800

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Filter-disc가 아닌 많은 구멍이 뚫려 있는 Perforated-plate
- 기본 Inner Joint Size : 24/40
- LK Labkorea



Plastic Buchner Funnel 플라스틱 부흐너 깔때기

Cat. No.	Model	용량	적용 필터	EA Price ₩
F07-59-111	148505	40 ml	Ø45 mm	24,500
F07-59-115	148510	70 ml	Ø55 mm	26,900
F07-59-119	148515	180 ml	Ø70 mm	27,500
F07-59-123	148520	285 ml	Ø80 mm	34,100
F07-59-127	148525	390 ml	Ø90 mm	48,200
F07-59-131	148530	810 ml	Ø110 mm	65,500
F07-59-135	148535	2.1 L	Ø160 mm	112,800
F07-59-139	148540	6 L	Ø240 mm	192,800

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Funnel의 상/하부 분리형으로 분해, 세척, 건조가 용이
- LK Labkorea



Optional Filter Disk

Cat. No.	Model	Mesh Size	적용 필터	EA Price ₩
F07-59-151	148605	1 mm	Ø45 mm	2,400
F07-59-155	148610		Ø55 mm	2,400
F07-59-159	148615		Ø70 mm	2,500
F07-59-163	148620		Ø80 mm	2,500
F07-59-167	148625		Ø90 mm	2,700
F07-59-171	148630		Ø110 mm	2,700
F07-59-175	148635		Ø160 mm	2,800
F07-59-179	148640		Ø240 mm	3,000

- 재질 : HDPE, 120°C 내열, Non-autoclavable
- 여과 효율 상승 및 여과지 손상 방지용

