

Adapter, for PFA or PP Compressing Fitting 피팅용 투입 어댑터, LukeGL®

- PFA Compression Fitting을 이용, Teflon 및 SUS Tubing 등과 연결 가능
- Borosilicate Glass 3.3
- LukeGL®

일자형, Straight ①

Cat. No.	Model	Cone	Connecting	EA Price ₩
A01-132-063	A1012C01	24/40	od1/2"	34,700
A01-132-067	A1012C02	29/42		35,900
A01-132-071	A1012C03	45/50		41,900
A01-132-075	A1013C01	24/40	od1/4"	34,700
A01-132-078	A1013C02	29/42		35,900
A01-132-081	A1013C03	45/50		41,900

곡형, 90° Bent ②

Cat. No.	Model	Cone	Connecting	EA Price ₩
A01-132-084	A1014C01	24/40	od1/2"	40,800
A01-132-087	A1014C02	29/42		42,000
A01-132-091	A1014C03	45/50		48,000
A01-132-095	A1015C01	24/40	od1/4"	40,800
A01-132-098	A1015C02	29/42		42,000
A01-132-101	A1015C03	45/50		48,000




Adapter, with female Luer Lock Luer Lock 팁 연결 어댑터

Cat. No.	Model	Joint	Type	EA Price ₩
A01-132-105	1055-01	14/20	Cone Type ①	129,400
A01-132-112	1055-04	24/40		141,000
A01-132-126	1055-10	14/20	Socket Type ②	135,800
A01-132-133	1055-13	24/40		147,100

- 재질 : PTFE, -200 ~ 280°C 내열, Autoclavable
- Adapter 가운데에 1/4" Female NPT Hole이 있음
- Compression Fitting이나 Pressure/Vacuum Gauge 연결용
- Chemglass

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
A01-132-656	1055-20	Female Luer Lock To Female Thread ③	43,100
A01-132-664	1164-N-02	Nut for od1/8" ④	5,800
A01-132-672	1164-F-02	Ferrule for od1/8" ⑤	3,400



Adapter, with Luer Outer Luer 팁용 어댑터

Cat. No.	Model	Cone	높이	EA Price ₩
A01-132-144	1055-A-14	14/20	50 mm	56,000
A01-132-148	1055-A-01	24/40	80 mm	64,600

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Chemglass

