

PTFE Buchner Funnel

테프론 부흐너 깔대기

Cat. No.	Model	용량	사용 필터	Size	EA Price ₩
F07-460-031	N1654-02	50ml	Ø 45mm	h94mm	548,600
F07-460-035	N1654-04	75ml	Ø 50mm	h117mm	537,500
F07-460-039	N1654-06	135ml	Ø 70mm	h142mm	667,200
F07-460-043	N1654-08	290ml	Ø 90mm	h165mm	1,047,200

- 재질 : PTFE , -200 °C to 250°C
- 용도 : Suspension of Solids
- 상/하 분리형 : Easy Cleaning, Easy Removal of filter Calce
- BOLA

- PTFE Filtering Disk는 424 Page 참조



Table Top Buchner Funnel, PE

탁상형 부흐너 깔때기

- 재질 : PE Body, HDPE Filter Plate
- Connection : id12.7 mm Tubing과 연결
- Filter Plate 일체형과 분리형으로 구분되어짐
- Filter Plate : 다공성 (Porous) 또는 Perforated
- 사용 온도 : up to 52°C
- Chemglass

with Porous Filter Plate

Cat. No.	Model	Plate Type	Funnel (id, mm)	전체 높이 (mm)	내부 높이 (mm)	Pore	EA Price ₩
F07-132-730	CG-1413-01	일체형	260	178	127	M	1,081,200
F07-132-734	CG-1413-02					C	1,066,300
F07-132-738	CG-1413-05		457	292	229	M	2,645,900
F07-132-742	CG-1413-06					C	2,645,900
F07-132-746	CG-1413-09					M	5,498,200
F07-132-750	CG-1413-10	610	330	267	C	5,498,200	
F07-132-754	CG-1413-03	분리형	260	178	127	M	1,081,200
F07-132-758	CG-1413-04					C	1,081,200
F07-132-762	CG-1413-07		457	292	229	M	2,645,900
F07-132-766	CG-1413-08					C	2,645,900
F07-132-770	CG-1413-11					M	5,498,200
F07-132-774	CG-1413-12		610	330	267	C	5,854,500

- Medium (M) : 45 ~ 90 µm, Coarse (C) : 90 ~ 130 µm
- Filter Plate 두께 : 6.4 mm

with Perforated Filter Plate

Cat. No.	Model	Plate Type	Funnel (id, mm)	전체 높이 (mm)	내부 높이 (mm)	EA Price ₩
F07-132-778	CG-1413-17	일체형	260	178	127	1,081,200
F07-132-782	CG-1413-19		457	292	229	2,286,000
F07-132-786	CG-1413-21		610	330	267	3,625,900
F07-132-790	CG-1413-18	분리형	260	178	127	1,081,200
F07-132-794	CG-1413-20		457	292	229	2,286,000
F07-132-798	CG-1413-22		610	330	267	3,625,900

- Perforation : Ø4.8 mm (at Ø11 mm)
- Filter Plate 두께 : 4.8 mm

Filter Paper

Cat. No.	Model	사용 Funnel id	Pore	PK Price ₩
F07-132-802	CG-1413-30	260 mm	3 µm	100ea 288,900
F07-132-806	CG-1413-31	457 mm		100ea 538,700
F07-132-810	CG-1413-32	610 mm		100ea 686,600

- Flow Rate : 90 ml / min
- Flow Rate : 90 ml / min 재질 : Creped White Cellulose
- 규격 : 3 µm pore / 0.3t

