

Porcelain Buchner Funnel


자제 부흐너 깔때기

Cat. No.	Model	용량	적용 필터	EA Price ₩
F07-128-111	237/1	35 ml	Ø45 mm	19,600
F07-128-115	237/2	70 ml	Ø55 mm	23,100
F07-128-119	237/3	120 ml	Ø70 mm	27,500
F07-128-123	237/4	240 ml	Ø90 mm	39,400
F07-128-127	237/5	400 ml	Ø110 mm	48,200
F07-128-131	237/6	600 ml	Ø125 mm	62,800
F07-128-135	237/7	1 L	Ø150 mm	82,800
F07-128-139	237/8	2 L	Ø185 mm	114,100
F07-128-143	237/9	4 L	Ø240 mm	196,300
F07-128-147	237/9a	7 L	Ø295 mm	556,100
F07-128-151	237/10	10 L	Ø325 mm	695,300

- 자기 재질로 내화학성 및 내열성이 뛰어남
- Max. Temp : 1,000°C, Density : 2.4 g/cm³
- Jipo / LabCeramic®

상세 규격

Model	D1	D2	H1	H2	Porosity
237/1	48 mm	38 mm	24 mm	43 mm	1 mm
237/2	62 mm	45 mm	30 mm	64 mm	
237/3	77 mm	58 mm	35 mm	64 mm	
237/4	97 mm	70 mm	40 mm	71 mm	
237/5	116 mm	95 mm	49 mm	83 mm	2 mm
237/6	130 mm	110 mm	52 mm	85 mm	
237/7	156 mm	130 mm	59 mm	96 mm	
237/8	192 mm	160 mm	74 mm	117 mm	
237/9	248 mm	220 mm	90 mm	117 mm	
237/9a	296 mm	270 mm	105 mm	130 mm	
237/10	334 mm	300 mm	120 mm	133 mm	



- Filter Adapter는 493 Page 참조



Stainless Steel Buchner Funnel

스테인레스 부흐너 깔대기

Cat. No.	Model	적용 필터	Top od	Stem od	EA Price ₩
F07-695-020	92-0571	Ø110 mm	Ø117 mm	Ø19 mm	526,100
F07-695-023	92-0572	Ø150 mm	Ø157 mm		704,300
F07-695-026	92-0573	Ø180 mm	Ø192 mm	Ø22 mm	755,300
F07-695-029	92-0574	Ø240 mm	Ø247 mm	Ø40 mm	1,205,100
F07-695-032	92-0575	Ø300 mm	Ø307 mm		1,374,800
F07-695-035	92-0576	Ø360 mm	Ø367 mm		1,553,000

- 재질 : Stainless Steel 304
- Sansyo



- Membrane Filter는 421 ~ 424 Page 참조

