

## Stainless Steel Funnel 스테인레스 깔때기

• 재질 : Stainless Steel 304 • Universal

### Funnel

Cat. No.	Model	Size	Stem Size	EA Price ₩
F07-215-101	904-11	∅30×h29 mm	∅11 mm	22,800
F07-215-105	904-12	∅45×h35 mm		26,300
F07-215-106	92-0560	∅87×h72 mm		32,200
F07-215-107	92-0559	∅115×h100 mm	∅15 mm	36,100
F07-215-108	92-0590	∅145×h110 mm		43,500

### Powder Funnel

Cat. No.	Model	Size	Stem Size	EA Price ₩
F07-215-109	904-15	∅148×h110 mm	∅27 mm	78,900
F07-215-113	904-16	∅207×h175 mm	∅48 mm	125,500



## Stainless Steel Funnel with Long Handle 손잡이형 스테인레스 깔때기

Cat. No.	Model	Dia x total length	Handle Size	EA Price ₩
F07-215-205	92-1555	∅48×L82 mm	45 mm	33,100
F07-215-210	92-1556	∅60×L90 mm	60 mm	34,000
F07-215-215	92-1557	∅76×L155 mm	70 mm	34,800
F07-215-220	92-1558	∅93×L220 mm	90 mm	36,500
F07-215-225	92-1559	∅120×L248 mm	120 mm	78,100

• 재질 : Stainless Steel • 재질 : Wire Type의 Handle 부착  
• Sansyo



92-1555



92-1557

## Stainless Funnel with Handle 손잡이형 스테인레스 깔때기

Cat. No.	Model	Overall Size	Stem Size	EA Price ₩
F07-102-210	4431	∅100×h120 mm	∅14.5×h60 mm	49,000
F07-102-220	4432	∅120×h140 mm	∅14.5×h60 mm	52,800
F07-102-230	4433	∅150×h170 mm	∅25×h90 mm	64,900
F07-102-240	80785200	∅200×h220 mm	∅20×h90 mm	55,800
F07-102-250	80785250	∅250×h260 mm	∅20×h90 mm	75,500

• 재질 : Stainless Steel 304 • LK Labkorea



## Stainless Steel Funnel 스테인레스 깔때기, Round-type

Cat. No.	Model	Size	Stem Size	EA Price ₩
F07-21-111	CYS5603	∅150×h130 mm	∅16×L32 mm	22,900
F07-21-118	CYS5604	∅200×h160 mm	∅32×L50 mm	33,700
F07-21-125	CYS5605	∅250×h170 mm	∅38×L52 mm	45,800

• 재질 : Stainless Steel 304 • LK Labkorea



## Stainless Steel Funnel for Filter Disc 스테인레스 깔때기, Mesh 교체형

Cat. No.	Model	Size	Stem Size	EA Price ₩
F07-695-010	92-1496	∅123×h87 mm	∅27 mm	26,900

• 재질 : Stainless Steel 304  
• 2종류의 Mesh Filter Disk와 같이 사용 가능  
• Mesh Filter Disk는 별도 구매해야함  
• Sansyo

### Optional Mesh Filter Disk

Cat. No.	Model	Mesh Size	Size	EA Price ₩
F07-695-013	92-1497	30	od98 mm	6,600
F07-695-016	92-1498	40		6,600

