

Separatory Funnel, Conical Type, Lenz

표준형 분액깔때기

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- PTFE Stopcock 및 PP Stopper 기본 포함
- Lenz

Without Graduation, 표준형 일반

Cat. No.	Model	용량	Stopper	Bore	EA Price ₩
F07-83-672	4004128	50 ml	19/26	2.5 mm	62,200
F07-83-678	4004137	100 ml			66,100
F07-83-684	4004149	250 ml	29/32	4.0 mm	79,600
F07-83-690	4004158	500 ml			95,100
F07-83-696	4004170	1 L		6.0 mm	127,200
F07-83-702	4004173	2 L			178,700

With Graduation, 표준형 눈금

Cat. No.	Model	용량	Stopper	Bore	눈금	EA Price ₩
F07-83-708	4006128	50 ml	19/26	2.5 mm	2 ml	62,800
F07-83-714	4006137	100 ml			5 ml	70,200
F07-83-720	4006149	250 ml	29/32	4.0 mm	10 ml	84,800
F07-83-726	4006158	500 ml			20 ml	103,100
F07-83-732	4006170	1 L		6.0 mm	50 ml	136,300
F07-83-738	4006173	2 L			192,800	



Separatory Funnel, Pear Type, LukeGL®

피어형 분액 깔때기


- Borosilicate Glass 3.3
- PTFE Stopcock 및 PP Stopper 기본 포함
- LukeGL®

Without Graduation, 피어형 일반

Cat. No.	Model	용량	Stopper	Bore	EA Price ₩
F07-83-201	4003128	50 ml	19/26	2.5 mm	41,000
F07-83-207	4003137	100 ml			42,000
F07-83-213	4003149	250 ml	29/32	4.0 mm	46,000
F07-83-219	4003158	500 ml			52,000
F07-83-225	4003170	1 L		6.0 mm	74,000
F07-83-231	4003173	2 L			109,000

With Graduation, 피어형 눈금

Cat. No.	Model	용량	Stopper	Bore	눈금	EA Price ₩
F07-83-237	4005128	50 ml	19/26	2.5 mm	2 ml	47,200
F07-83-243	4005137	100 ml			5 ml	48,300
F07-83-249	4005149	250 ml	29/32	4.0 mm	10 ml	52,900
F07-83-255	4005158	500 ml			59,800	
F07-83-261	4005170	1 L		6.0 mm	20 ml	85,100
F07-83-267	4005173	2 L			125,400	



Separatory Funnel, Squibb Type, LukeGL®

장형 분액깔때기

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- PTFE Stopcock 및 PP Stopper 기본 포함
- LukeGL®

Without Graduation, 장형 일반

Cat. No.	Model	용량	Stopper	Bore	EA Price ₩
F07-83-744	4014128	50 ml	19/26	2.5 mm	34,200
F07-83-750	4014137	100 ml			35,100
F07-83-756	4014149	250 ml	29/32	4.0 mm	38,000
F07-83-762	4014158	500 ml			42,300
F07-83-768	4014170	1 L		6.0 mm	61,300
F07-83-774	4014173	2 L			90,600

With Graduation, 장형 눈금

Cat. No.	Model	용량	Stopper	Bore	눈금	EA Price ₩
F03-83-780	4016128	50 ml	19/26	2.5 mm	2 ml	37,300
F03-83-786	4016137	100 ml			5 ml	38,700
F03-83-792	4016149	250 ml	29/32	4.0 mm	10 ml	42,200
F03-83-798	4016158	500 ml			47,000	
F03-83-804	4016170	1 L		6.0 mm	20 ml	65,500
F03-83-810	4016173	2 L			95,800	

