

Amber Glass Separatory Funnel

갈색 분액 깔때기


- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- PTFE Stopcock 및 PE Stopper 기본 포함
- Lenz

Conical Type ①

Cat. No.	Model	용량	Stopper	Bore	EA Price ₩
F07-83-780	4008128	50 ml	19/26	2.5 mm	135,500
F07-83-782	4008137	100 ml			151,100
F07-83-784	4008149	250 ml	29/32	4.0 mm	190,800
F07-83-786	4008158	500 ml			221,200
F07-83-788	4008170	1 L		6.0 mm	287,800
F07-83-790	4008173	2 L			402,700

Squibb Type ②

Cat. No.	Model	용량	Stopper	Bore	EA Price ₩
F07-83-792	4018128	50 ml	19/26	2.5 mm	141,900
F07-83-794	4018137	100 ml			144,200
F07-83-796	4018149	250 ml	29/32	4.0 mm	172,800
F07-83-798	4018158	500 ml			207,700
F07-83-800	4018170	1 L		6.0 mm	256,900
F07-83-802	4018173	2 L			300,000



Large Separatory Funnel, Round Type

대용량 분액깔때기, LukeGL®

Cat. No.	Model	용량	Stopper	Bore	EA Price ₩
F07-83-816	L4020128	3 L	34/45	6 mm	203,100
F07-83-819	L4020133	4 L			241,000
F07-83-822	L4020137	5 L			361,400
F07-83-825	L4020142	6 L	45/50	8 mm	483,100
F07-83-828	L4020149	10 L			680,200
F07-83-834	L4020158	20 L			860,600

- Borosilicate Glass 3.3
- PTFE Stopcock 및 PP Stopper 기본 포함
- LukeGL®



Plastic Separatory Funnel with Screwcap

플라스틱 분액 깔때기, 스크류 캡 포함

Cat. No.	Model	용량	EA Price ₩
F07-339-501	27159-0125	125 ml	28,400
F07-339-503	27159-0250	250 ml	40,100
F07-339-505	27159-0500	500 ml	46,100
F07-339-507	27159-1L	1 L	51,800

- 재질 : PP Body & Screwcap, -10 ~ 135°C 내열, Autoclavable
- PTFE Stopcock, Mold 눈금
- LabPlasti®



Separatory Funnel with Stopper

분액 깔때기, Stopper 포함

Cat. No.	Model	용량	EA Price ₩
F07-339-511	27157-0125	125 ml	21,050
F07-339-517	27157-0250	250 ml	33,500
F07-339-523	27157-0500	500 ml	38,100
F07-339-525	27157-1L	1 L	42,700

- Body 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- PTFE Stopcock, PP Stopper, Mold 눈금
- Stopcock와 몸체를 분리 할 수 있어 세척이 용이
- LabPlasti®

