

Teflon Separatory Funnel 테프론 분액 깔때기

Cat. No.	Model	Capacity	EA Price ₩
F07-23-214	TSF-100	100 ml	223,100
F07-23-232	TSF-200	200 ml	334,800
F07-23-226	TSF-500	500 ml	569,300

- 재질 : PFA Body (FEP 500 ml Only), PTFE Cock
- Universal



FEP Dropping Funnel with Joint FEP 테프론 드로핑 깔때기, 하부 조인트

Cat. No.	Model	용량	Joint	EA Price ₩
F03-460-011	B285-01	125 ml	29/32	1,021,000
F03-460-015	B285-02	250 ml		1,162,500
F03-460-019	B285-03	500 ml		1,512,300

- 재질 : FEP, 0 °C to 110°C 내열
- PTFE Stopcock with 29/32 Cone
- Ground Joint의 Easy Locking과 Unlocking을 위한 Safety Nut 포함
- BOLA



Separatory Funnel Rack은 840 Page 참조

