

Program Muffle Furnace 프로그램식 전기로, 1100°C & 1200°C

배송비 별도

Cat. No.	Model	Capacity	Heater	Max. Temp.	EA Price ₩
F08-101-010	KL0311HT	3 L	1.2 kW	1100°C	4,563,700
F08-101-020	KL0511HT	5 L	2.4 kW		5,799,300
F08-101-030	KL0911HT	10 L	3.0 kW		6,058,300
F08-101-040	KL1511HT	15 L	3.5 kW		7,294,000
F08-101-050	KL0312HT	3L	1.2 kW	1200°C	5,560,100
F08-101-060	KL0512HT	5 L	2.4 kW		6,696,000
F08-101-070	KL0912HT	10 L	3.0 kW		7,772,200
F08-101-080	KL1512HT	15 L	3.5 kW		8,071,200

- 일반 실험실에서 널리 사용 되는 Furnace임
- Stainless Steel Case로 내구성이 높음
- 고 성능의 Heating Wire는 Ceramic Plate에 몰딩되어 있어 물리적 충격에 강함
- Gas 배기구는 Furnace 뒷면에 있음
- Eurotherm 3208 Program Controller (최대 8개까지의 Segment)
- Adjustable Temperature Limiter (Option)
- 약세사리는 약세사리란 참조
- Thermconcept



Cat. No.	F08-101-010	F08-101-020	F08-101-030	F08-101-040
Model	KL0311HT	KL0511HT	KL0911HT	KL1511HT
Temp. Range	Max. 1100°C			
Cat. No.	F08-101-050	F08-101-060	F08-101-070	F08-101-080
Model	KL0312HT	KL0512HT	KL0912HT	KL1512HT
Temp. Range	Max. 1200°C			
Temp. Controller	Eurotherm 3208 Program Controller			
용량	3 L	5 L	10 L	15 L
Heater 용량	1.2 kW	2.4 kW	3.0 kW	3.5 kW
Size (w×d×h, mm)	Internal	190×140×100	240×170×130	240×260×170
	External	440×445×490	460×465×520	495×540×560
Weight	20 kg	35 kg	45 kg	50 kg
Power	230V, 1/N			

Compact Muffle Furnace, 1000°C 소형 전기로

Cat. No.	Model	용량	Heater	Max. Temp.	EA Price ₩
F08-101-092	KLS0310	3 L	1.5 kW	1000°C	3,945,900

- User-friendly Swing door
- 벽면 전체에 열선이 분포 되어 있어 Heating Cycle이 빠름
- Gas 배기구는 Furnace 뒷면에 있음
- PC442T Controller
- Microprocessor controllers ensure precision furnace management for both simple and complex processes
- 약세사리는 약세사리란 참조
- Thermconcept

Cat. No.	F08-101-092	
Model	KLS0310	
Temp. Range	Max. 1000°C	
Temp. Controller	PC442T Controller	
용량	3 L	
Heater 용량	1.5 kW	
Size (w×d×h, mm)	Internal	140×200×110
	External	320×350×410
Weight	20 kg	
Power	230V, 1/N	

