

## Laboratory Chamber Furnace, Program 프로그램식 전기로

배송비 별도

Cat. No.	Model	Capacity	Heater	Max. Temp.	EA Price ₩
F08-101-188	KLS0512	7 L	2.0 kW	1200°C	6,676,100
F08-101-191	KLS1012	10 L	3.6 kW		7,772,200
F08-101-194	KLS1512	15 L			7,971,600
F08-101-197	KLS3012	30 L	6.0 kW		10,940,800
F08-101-202	KLS4512	45 L	7.5 kW		11,339,500
F08-101-205	KLS0513	7 L	2.0 kW	1300°C	7,094,700
F08-101-208	KLS1013	10 L	3.6 kW		8,091,100
F08-101-211	KLS1513	15 L			8,708,900
F08-101-214	KLS3013	30 L	6.0 kW		11,837,700
F08-101-217	KLS4513	45 L	7.5 kW		14,009,900

- Heating on Both Side : 양면 가열식
- Insulation Made of High Grade Ceramic Fibre with Low Thermal Mass
- 상/하 개폐식 Door • Gas 배기구는 Furnace 뒷면에 있음
- Double-Walled Housing : Low Outer-casing Temperature
- Ceramic Bottom Plate (Included)
- PC442T Controller
  - Microprocessor Controllers Ensure Precision Furnace Management for Both Simple and Complex Processes
- Adjustable Temperature Limiter (Option)
- Thermconcept



### Max. Temp 1200°C

Cat. No.	F08-101-188	F08-101-191	F08-101-194	F08-101-197	F08-101-202	
Model	KLS0512	KLS1012	KLS1512	KLA3012	KLS4512	
Temp. Controller	PC442T Controller					
용량	7 L	10 L	15 L	30 L	45 L	
Size (w×d×h)	Internal	200×250×140 mm	200×250×200 mm	220×300×230 mm	280×380×280 mm	300×500×300 mm
	External	550×580×650 mm	560×640×720 mm	560×640×720 mm	590×690×790 mm	660×810×820 mm
Heater 용량	2.0 kW	3.6 kW		6.0 kW	7.5 kW	
Power	230V, 1/N			400V, 3/N		

### Max. Temp 1300°C

Cat. No.	F08-101-205	F08-101-208	F08-101-211	F08-101-214	F08-101-217	
Model	KLS0513	KLS1013	KLS1513	KLS3013	KLS4513	
Temp. Controller	PC442T Controller					
용량	7 L	10 L	15 L	30 L	45 L	
Size (w×d×h)	Internal	200×250×140 mm	200×250×200 mm	220×300×230 mm	280×380×280 mm	200×500×300 mm
	External	550×580×650 mm	560×640×720 mm	560×640×720 mm	590×690×790 mm	660×810×820 mm
Heater 용량	2.0 kW	3.6 kW		6.0 kW	7.5 kW	
Power	230V, 1/N			400V, 3/N		

## Accessory for "Thermconcept" Furnace "Thermoconcept" 전기로용 악세서리

### Exhaust with Fan

Cat. No.	Model	Description	EA Price ₩
F08-101-950	KL-EF	KL-Type 전기로용 Vent Fan	1,016,400



### High Temp. Glove / 고온용 장갑

Cat. No.	Model	Description	EA Price ₩
F08-101-980	41103	Max. 600°C	148,700
F08-101-990	41101	Max. 900°C	238,300

