

## 1600°C High Temp. Program Furnace 프로그램식 고온 전기로

배송비 별도

Cat. No.	Model	Capacity	Heater	Max. Temp.	EA Price ₩
F08-101-755	KLC0514	5 L	3.8 kW	1400°C	15,751,400
F08-101-758	KLC1014	9 L	5.0 kW		22,534,800
F08-101-761	KLC1514	15 L	7.0 kW		30,697,900
F08-101-764	KLC3014	30 L	11.5 kW		43,460,000
F08-101-767	KLC0515	5 L	4.8 kW	1500°C	17,521,900
F08-101-773	KLC1015	9 L	6.6 kW		24,144,500
F08-101-776	KLC1515	15 L	7.5 kW		37,251,400
F08-101-782	KLC0516	5 L	5.2 kW	1600°C	19,499,500
F08-101-785	KLC1016	9 L	6.9 kW		28,743,400
F08-101-788	KLC1516	15 L	8.0 kW		36,101,700
F08-101-791	KLC3016	30 L	13.5 kW		54,957,300

- Heating on Both Side : 2면 가열식 (Chamber내 균일한 온도 분포를 나타냄)
- Powerful SiC-rod Heating Element
- Double-walled Housing with Rear-ventiation to Ensure Low Outer Casing Temperatures
- Ceramic Bottom Plate (Included)
- Gas 배기구는 Furnace 윗면에 있음
- Paralled-guided Door Moving Up-wards
- Eurotherm 3208 Program Controller
- Programmable PID Controller with Realtime Clock
- Simple, Straightforward Operation
- 1 Program - 8 Segments (4 Ramp and 4 Holding Time)
- Optional RS232 or RS485 interface for PC connection (iTools 소프트웨어)
- Adjustable Temperature Limiter (Option)
- 약세사리는 515 Page 참조
- **Thermconcept**



### Technical Data

Cat. No.	F08-101-755	F08-101-758	F08-101-761	F08-101-764
Model	KLC0514	KLC1014	KLC1514	KLC3014
Temp. Range	Max. 1400°C			
용량	5 L	9 L	15 L	30 L
Heater 용량	3.8 kW	5.0 kW	7.0 kW	11.5 kW
Size	Internal	w150×d250×h140 mm	w200×d250×h180 mm	w220×d310×h220 mm
	External	w550×d580×h650 mm	w560×d670×h720 mm	w590×d690×h790 mm
Power	400V, 3/N			

Cat. No.	F08-101-767	F08-101-773	F08-101-776
Model	KLC0515	KLC1015	KLC1515
Temp. Range	Max. 1500°C		
용량	5 L	9 L	15 L
Heater 용량	4.8 kW	6.6 kW	7.5 kW
Size	Internal	w150×d250×h140 mm	w200×d250×h180 mm
	External	w550×d580×h650 mm	w560×d640×h720 mm
Power	400V, 3/N		

Cat. No.	F08-101-782	F08-101-785	F08-101-788	F08-101-791
Model	KLC0516	KLC1016	KLC1516	KLC3016
Temp. Range	Max. 1600°C			
용량	5 L	9 L	15 L	30 L
Heater 용량	5.2 kW	6.9 kW	8.0 kW	13.5 kW
Size	Internal	w150×d250×h140 mm	w200×d250×h180 mm	w220×d310×h220 mm
	External	w550×d580×h650 mm	w560×d670×h720 mm	w590×d690×h790 mm
Power	400V, 3/N			
Temp. Controller (공통)	Eurotherm 3208 Program Controller			