

Ductless Storage Cabinet 필터식 시약장, Chemical 필터 포함

배송비 별도

Cat. No.	Model	Size (w×d×h, mm)	Door	EA Price ₩
F09-138-100	NS800S	w800×d525×h2362mm	2 door	4,626,600
F09-138-140	NS1600	w1600×d525×h2362mm	4 door	9,181,300

- 필터를 통해 시약장 안에 보관된 화학물질들의 해로운 냄새 제거 및 24시간 실험실의 공기정화
- Organic Chemical용 필터 기본 장착
- 요청에 의해 필터 타입 변경 가능 (Chemical 종류에 따른 다양한 필터 : Organic chemical / Acid / Amonia / Formaldehyde)
- 디스플레이 창을 통해 현재 상태 안내 : Fan 스피드 / 필터 배출 농도 / 내부 공기 상태 / 내부 온도
- 필터에서 배출 되는 공기의 Chemical 의 농도를 실시간으로 보여 줌으로써 필터 교체 여부를 알 수 있음 (검출 한계는 0.1ppm 이지만 환경에 대한 고객의 다양한 요구 사항에 따라 경보 값을 설정할 수 있음)
- Fan 스피드 조절 가능 : 최상의 흡착 효과와 안전성을 보장하려면 팬의 속도를 적절하게 조절 해야함
- 600여 종 이상의 다양한 화학물질 필터링
- Hua Feng

Cat. No.	F09-138-100	F09-138-140
Model	NS800S	NS1600
Volume	160 L	320 L
	1L 시약병 약130개	1L 시약병 약260개
Material	Exterial	내산성 Epoxy Coated Steel 1.2t
	Shelf	Polypropylene Plastic Tray
	Door	High Quality Acrylic 5t
Size	w800×d525×h2362mm	w1600×d525×h2362mm
Adjustable Shelf	8	16
	Max. Load : 각 50Kg	
공기처리량	240 m3/h	
Noise	< 42	< 46
Power	220V, 50/60Hz	

Filter for Ductless Storage Cabinet

Cat. No.	Model	Description	EA Price ₩
F09-138-710	HF4YJ	Filter Organic Chemical	502,200
F09-138-720	HF4TY	Filter Acid Chemical	542,800
F09-138-730	HF4FM	Filter Formaldehyde	602,500
F09-138-740	HF4AM	Filter Amonia	578,600
F09-138-750	HF4HP	HEPA Filter	542,800

