

Natural Rubber Balloon 고무 풍선

Cat. No.	Model	용량	두께	PK Price ₩
G01-161-070	52040	7 L	12 mil (0.03 t)	100ea 31,900

- 재질 : 100% Pure Natural Rubber
- 100 ea/pk
- 색상 : 랜덤색상
- Betallatex



Gas Balloon Flask 가스 저장 플라스크

Cat. No.	Model	용량	Needle / Valve	EA Price ₩
G01-132-103	1800-01	100 ml	0 ~ 4 mm Pass	382,500
G01-132-109	1800-02	250 ml		406,600
G01-132-115	1800-03	500 ml		414,300

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 0 ~ 4 mm High Vacuum Valve
- Chemglass



Gas Sampling Tube 가스 샘플링 튜브

Cat. No.	Model	용량	Teflon Cock	EA Price ₩
G01-132-127	1808-01	125 ml	Bore 2 mm	344,000
G01-132-133	1808-02	250 ml	Bore 4 mm	390,600


- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Sampling Pore : od8 mm용 Sleeve septa 사용
- Chemglass



PP Gas Sampling Bulb 가스 샘플링 벌브

Cat. No.	Model	용량	Bulb Size	EA Price ₩
G01-46-031	021010	250 ml	Ø51×L179 mm	69,400

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- 2-way Teflon Cock
- LabPlasti®



Syphon 사이폰

Cat. No.	Model	Size	Tube id	EA Price ₩
G01-46-333	660010	Ø60×L210 mm	6.0 mm	13,000
G01-46-337	660020	Ø80×L280 mm	12.5 mm	31,000

- 재질 : LDPE, -50 ~ 80°C 내열
- LabPlasti®



Gas leakage Checker 가스 누출 체크액

Cat. No.	Model	용량	방식	EA Price ₩
G10-64-615	F7015-42	8 oz	Liquid Type	15,500

- 냉각제, 프로판가스, 도시가스 등의 누출시 배관이나 배관 접속부에 사용함
- Flon Industry

